

HEALTH & BEAUTY

COMPLETE CATALOGUE EN 1/2026

CONTENTS

Nail care

X-Series	6
Files for professional nail design	8
Files for home use / travel	16
Ultrasonic cleaning	20
UV lamps + dust extraction	21

Health check

Air purifiers	27
Oral hygiene	29
Blood pressure monitors	30
Thermometers	34
Oximeters	35
Blood glucose monitors	37

Therapy

TENS/EMS	38
----------------	----

Wellness

Massage	40
Foot spa	41
Intimate care	42
Heat	43
Room air / aroma diffuser	50
Aromatic oils.....	54

Shoes, soles and accessories

Introduction to Promed comfort shoes	56
Overview of Promed shoes and indications.....	60
Promed size tables.....	62
Promed bandage shoes	63
Huber soles by Promed	72

Accessories

Grinding bits, handpieces, pedals etc.	74
Electrodes etc.	92
Cuffs, mouthpieces, masks, etc.	94
Oral hygiene nozzles, blood sugar accessories	95

Glossary	96
----------------	----

Advertising materials	98
-----------------------------	----

Legal notice	99
--------------------	----



Advantages of Promed files

We offer the latest technology and the best service from one single source.

High quality for durable products.

Our highly qualified staff are at your service if you have any questions.

● SINCE 1982
● **GERMAN**
● **BRAND**

Professional appliances for skin and nail care

As a leading manufacturer of technically sophisticated and tried-and-tested professional milling tools, Promed develops, researches and designs its appliances in close collaboration with both the users and the medical technology engineers.

We translate our more than 40 years of experience into useful developments for you that make the work you do for your customers much easier and quicker and more comfortable.

Promed was the first manufacturer to develop technical innovations such as ASC (automatic speed control) and the tried-and-tested OPC (overload protection control).

Promed appliances are subjected to both regular tests by recognised testing institutes as well as the highest German and global safety standards. In addition to outstanding technical innovations, our strengths also lie in our unique after sales service. If you wish, we can offer you a loaned appliance (only in Germany) for the annually-recommended customer service of your appliance.

We are represented at the following beauty trade fairs in Germany:



promed

Technology that makes the difference

APE – automatic power electronics

more power for effective work. The intelligent electronics provide the motor with additional power, also at low speeds.

ASC – automatic speed control

The digital automatic speed control notices the loss of speed as soon as pressure is applied to the drill bit during the filing and automatically compensates for it. The appliance maintains its set speed. The full power is always available to you. This allows you to work especially swiftly.

OPC – overload protection control

The unique overload protection automatically detects any overloading of the respective motor and therefore prevents the appliance from being damaged. This is available for handpieces and in some instances for suction turbines.

Overload protection for the handpiece motor:

If the handpiece is blocked or overloaded, the handpiece and/or the suction turbine will switch off.

Overload protection for the suction turbine:

If the suction turbine does not receive sufficient cooling and the air temperature in the suction turbine casing exceeds a critical value, the suction turbine and the handpiece will switch off.

Safety stop

If excessive speed and pressure are applied during the sanding, the electronics will detect this and will switch off the appliance. In this case, you will need to set the appliance to zero (centre position) with the rotation direction switch/controller in order to resume use of the appliance.

Surfaces

Nail polish remover-resistant membrane in an elegant silver metallic look.

Performance

A higher performance for more torque is achieved through the intelligent, microprocessor-controlled electronics.

All-in-one controller

Large, practical and easy to operate, it also allows you to choose the speed and direction of rotation.

Handpieces

New development made from an aluminium/stainless steel combination in a lightweight design with precision ball bearings. Especially high-performance micro-motors. Exceptionally ergonomic and vibration-free in operation. Non-metal handpieces are made of special plastic.

Promed – benefit from our know-how.



Hand care



Foot care



Artificial nails



ASC (automatic speed control)



OPC (overload protection control)



APE (automatic power electronics)



Safety stop



Cold light/LED light



Battery operated

NEW!



Nailfan Mini-X



UVL 60-X



The File 3030-X

X-Series: Wireless. Powerful. Durable.

The new battery-powered device family – designed for everyday studio use.



Nailfan Mini-X Professional suction device with battery

- **Brushless motor** Durable, quiet and low maintenance
- **10 suction levels** Can be finely adjusted to suit treatment and environment
- **Strong suction performance** Clean workplace
- **Up to 8 hours runtime** Mobile, ready again quickly

Details and technical data on page 22



UVL 60-X High-power LED lamp with battery

- **42 extra-large high-performance LEDs** For even results
- **Gold-coated LED contacts** Optimised conductivity, long service life
- **Top & rear operation** Perfect for desk in studio setup
- **Timer** 10/30/60/90 seconds
- **High-Power/Low-Heat** Flexible in use

Details and technical data on page 23



The File 3030-X Professional battery-powered file

- **Metal casing** Stable and high-quality for everyday studio use
- **Brushless motor** Efficient, durable, consistent performance
- **Up to 30,000 rpm** Precise, infinitely adjustable
- **Low-vibration handpiece** Ergonomic, low fatigue

Details and technical data on page 13

All X-Series devices are lithium battery powered. They are designed for the most demanding studio requirements, offering high mobility, ease of use and long-lasting components.

the file - 4040

State-of-the-art file with suction
with brushless handpiece and suction motor



INNOVATION AND TECHNOLOGY

"Made in Germany"

- + **1 Consistently high suction capacity!**
- + **2 Very quiet!**
- + **3 Lightweight, ergonomic handpiece!**



- + **1 Dual suction motor (brushless):**
The file 4040 is equipped with two brushless suction turbines. The high vacuum created by these ensures a consistently high suction capacity. During the milling, the extracted fine particulate goes from the handpiece head, through the hose and into the dust bag more quickly and easily. The brushless motors mean that there is no carbon dust soiling.

- + **2 Especially quiet:**
Despite the high suction performance, the special noise insulation allows for pleasant working conditions.

- + **3 Ergonomic, lightweight handpiece (approx. 100 g) with automatic collet, suction hose and (filter) drawer:**

Very lightweight, ergonomically shaped plastic handpiece with high-performance motor. High-quality finish. Ideal for working on finger and toe nails without damaging your joints. The micro-motor in the handpiece has a speed of up to 40,000 rpm. The standby button ensures you can quickly switch the handpiece on and off while working. Powerful 360-degrees suction directly at the milling tool.

The automatic collect allows for the quick and easy changing of bits. Easy cleaning/disinfection of the handpiece and the highly flexible suction hose.

- + **4 High-tech colour touch display (4.3 inches) with brushed front membrane:**

Large, coloured, easy-to-use touch display. The handpiece speed can be adjusted up to max. 40,000 rpm. Suction power adjustable between 6 levels (Eco to Max). 3 memory buttons for repeating work steps. Display of direction of rotation, speed, suction capacity, safety information as well as service and filter replacement.

Suction adjustable at 6 levels
(Eco - Max)

Up to 40,000 rpm
(with -/+ buttons or
touch/slide
adjusters)



Clockwise/anticlockwise rotation

Memory buttons



+ CLEAN AND CONSCIENTIOUS WORK

To protect your health and your environment, the fine particulate is extracted directly at the handpiece and goes into the filter bag during work.

The amount of dust in the studio is reduced to the minimum and affects neither your environment nor your health.



- + **5 Dust monitoring with safety switch-off:**
The built-in dust sensor monitors the accumulation of dust in the appliance. In case of dust in the appliance – e.g. caused by a defective dust bag, the appliance switches off automatically to prevent damage to the motor.

- + **6 Overheating protection (motor + handpiece):**
The file is equipped with an automatic switch-off in case of overheating of the suction or handpiece motor. If the appliance or the handpiece gets too hot, the motors switch off automatically.

- + **7 Expelled air filter to prevent fine particulate:**
To protect the environment and your health, the integrated extracted air filter cleans the expelled air of fine particulate.



Dynamic, flat design
– for your studio



Expelled air grille

- Large display (4.3 inches) with touch and slider function, digital colour display
- Infinitely adjustable up to 40,000 rpm (clockwise/anticlockwise rotation)
- 3 memory buttons for individual work steps
- Indicator for filter bag replacement
- Overheating protection: temperature-based safety motor shut-down
- Strong suction performance – very quiet: adjustable with 6 levels (Eco - Max)
- Extremely lightweight, ergonomic handpiece (approx. 100 g) with standby button (dust-protected)

- Flexible suction hose (with a length of approx. 140 cm)
- Including hang-up bracket for easy working
- Maintenance indicator / maintenance interval indicator
- Including 2 filter bags (approx. 14x18 cm), exchangeable, tear-proof, extra large
- Optional: Foot pedal (on/off), see accessories, p. 90

TECHNICAL DATA
Dimensions: (W x H x D) 265 x 140 x 330 mm • Weight: approx. 4.5 kg • Mains voltage: 110 - 230 VAC, 50/60 Hz • Speed: up to 40,000 rpm • Size of the display: 4.3 inches • Length of mains cable: approx. 200 cm • Hose length: approx. 140 cm Volume: suction capacity Eco = approx. 40 dB(A); suction capacity max. = approx. 54 dB(A)

ACCESSORIES



Filter bags 1 x (white)
+ 1 x in appliance



Hose hang-up bracket

OPTIONAL ACCESSORIES



Foot switch (on/off)



Handpiece tray

Art. 204036 | EAN 40 43641 20052 7

the file - 4030-SX2

For your health!

Modern file with suction

**Very low dust exposure
due to direct suction at
the grinding bits. Perfect for
nail design and foot care.**



- Membrane keypad (soft keys) with digital display
- Infinitely adjustable up to 30,000 rpm (clockwise/anticlockwise rotation)
- 3 memory keys for individual speeds
- Indicator for filter bag replacement
- Temperature-based safety shut-down in the case of overload
- Strong suction performance (infinitely adjustable 350 watt turbine), very quiet
- Extremely lightweight, ergonomic handpiece (95 g) with standby switch

- Flexible suction hose including hang-up bracket
- Including 5+1 filter bags (14x18 cm), exchangeable, tear-proof, extra large
- Including handpiece holder
- Optional: Foot pedal (on/off), and electronic rocker pedal (see accessories, p.90)

TECHNICAL DATA
Dimensions: (W x H x D) 225 x 170 x 230 mm · Weight: approx. 2.6 kg · Mains voltage: 230 VAC/50 Hz, 350 VA (115 VAC/60 Hz, 350 VA optional) · Turbine: 350 watts · Noise: 62 dBA, speed: up to 30,000 rpm

Product video:



promed


Art. 205050 | EAN 40 43 641 20061 9

the file - 5040-SX2

High-performance professional spraying appliance
for foot care using the wet technique

With colour display!

**All types of milling work in the area of
footcare. Nail and callus treatment.**



- 4 memory keys for programming individually-set speeds
- High-performance micro-pump for the finest spray mist (adjustable)
- Automatic jet cleaning when switching on and off
- Automatic pump cleaning
- Handpiece max. 40,000 rpm (clockwise/anti-clockwise rotation)
- Ergonomic handpiece with integrated spraying appliance
- Digital display of the recommended grinding bit

- Overload protection for the handpiece motor
- Standby key for all functions
- Standby rocker switch on handpiece
- Water reservoir (410 ml) – easy to refill with an attachment system
- Optional: Foot pedal (on/off), and electronic rocker pedal (see accessories, p.90)

TECHNICAL DATA
Power supply: 230 VAC/115 VAC, 50/60 Hz, 105 VA · 207 x 163 x 245 mm (L x W x H) · Pump capacity: 550 mbar, 8 ml/min. · Weight: Handpiece 112 g, appliance overall approx. 3.2 kg

Vacuum performance and number of revolutions




View of the keypad and filter replacement display



✓ ALLOWS FOR HEALTH-CONSCIOUS WORK

To prevent risks to your health, fine particulate is removed directly at the handpiece as you work. The accumulation of dust in the studio is reduced to a minimum. The practical guiding of the suction hose is ensured by our unique mounting bracket on the appliance.

ACCESSORIES



OPTIONAL ACCESSORIES



Foot switch (on/off) or foot pedal (rocker pedal)

NEW FEATURES OF THE APPLIANCE (compared to predecessor model)

① Automatic jet cleaning:

When the appliance is switched on, it starts by blowing only air through the nozzles in order to remove any dirt. When the appliance is switched off, the spray pump continues running for several seconds in order to clean the nozzle and/or the hose system.

② Spray setting

The file is equipped with a high-performance micro-pump which generates very fine spray. It offers the possibility to individually adjust the spray pattern according to the liquid used.



Scale for setting the spray pattern



Art. 117510 | EAN 40 43641 110130

The File - 725

Professional file for studio and home –
precise file work with professional comfort

NEW!

For professional manicures
and pedicures



- Speed infinitely adjustable up to 25,000 rpm
- Ergonomic, powerful, vibration-free stainless steel handpiece (160 g), special design for foot care
- Professional rotary collet
- 10 LEDs for speed display

- Safety stop
- ASC (automatic speed control)
- All-in-one controller
- Rotary knob with brightly sparkling jewel (diamond shape)

TECHNICAL DATA

Speed: Infinitely adjustable up to 25,000 rpm • Power supply: 100 - 240 VAC, 50/60 Hz, 1300 mA; Sec. 24 VDC, 1.5 A • Weight: 270 g • Weight of handpiece: 160 g • Dimensions: 150 x 135 x 75 mm (D x W x H)

promed

Art. 201030 | EAN 40 43641 110147

The File - 825

Professional file for studio and home – precise file
work with professional comfort (with foot pedal)

NEW!

For professional manicures
and pedicures



- Speed infinitely adjustable up to 25,000 rpm
- Ergonomic, powerful, vibration-free stainless steel handpiece (160 g)
- Professional rotary collet
- 10 LEDs for speed display
- 1 LED for foot pedal display
- ASC (automatic speed control)

- Safety stop
- All-in-one controller
- Rotary knob with brightly sparkling jewel (diamond shape)
- Foot pedal (on/off) (see accessories p.90)

TECHNICAL DATA

Speed: Infinitely adjustable up to 25,000 rpm • Power supply: 100 - 240 VAC, 50/60 Hz, 1300 mA; Sec. 24 VDC, 1.5 A • Weight: 275 g • Weight of handpiece: 160 g • Dimensions: 150 x 135 x 75 mm (D x W x H)

promed

Art. 203030 | EAN 40 43641 200640

the file - 3030-X

Professional file running on battery power
Professional quality in an elegant metal casing

NEW!

- Metal casing with integrated bit & handpiece support – robust, elegant & practical
- Extremely quiet working – hardly any vibrations, maximum smoothness
- Brushless high-performance motor – state-of-the-art technology for precision & long service life
- Infinitely adjustable up to 30,000 rpm – full control for all applications
- Low-vibration stainless steel handpiece – ergonomic and low fatigue

TECHNICAL DATA

• Station dimensions: approx. 100 x 100 x 40 mm • Handpiece dimensions: Ø approx. 18 x 139.5 mm
• Battery capacity: 2600 mAh • Station weight: approx. 440 g • Handpiece weight: approx. 138 g
• Power supply: 5 V 2 A DC • Speed: 30,000 rpm • Power: 78 W • Precise collet system

Art. 202025 | EAN 40 43641 20053 4

the file - 2020

The latest generation
of files for continuous work



For professional manicures
and pedicures



- Speed infinitely adjustable up to 20,000 rpm
- Ergonomic, powerful, vibration-free stainless steel handpiece (160 g), special design for foot care
- Professional rotary collet
- 5 LEDs for speed display
- 1 LED for foot pedal display
- Safety stop



- ASC (automatic speed control)
- All-in-one controller
- Rotary knob with brightly sparkling jewel (diamond shape)
- Foot pedal optionally available (on/off) (see accessories, p. 90)

TECHNICAL DATA
Speed: Infinitely adjustable up to 20,000 rpm · Power supply: 100 - 240 VAC, 50/60 Hz, 800 mA; Sec. 24 VDC, 1.5 A · Weight: 450 g · Weight of handpiece: 160 g · Dimensions: 130 x 200 x 100 mm (D x W x H)

Art. 201030 | EAN 40 43641 20036 7

the file - 1030

The latest generation of files for
professional nail design



For the professional treatment
of artificial and natural nails and
nail design in continuous use



- Speed infinitely adjustable up to 30,000 rpm
- Ergonomic, powerful, vibration-free stainless steel handpiece (160 g)
- Professional rotary collet
- 5 LEDs for speed display
- 1 LED for foot pedal display
- ASC (automatic speed control)



- Safety stop
- All-in-one controller
- Rotary knob with brightly sparkling jewel (diamond shape)
- Foot pedal optionally available (on/off) (see accessories, p. 90)

TECHNICAL DATA
Speed: Infinitely adjustable up to 30,000 rpm · Power supply: 100 - 240 VAC, 50/60 Hz, 800 mA; Sec. 24 VDC, 1.5 A · Weight: 450 g · Weight of handpiece: 160 g · Dimensions: 130 x 200 x 100 mm (D x W x H)

Art. 116510 | EAN 40 43641 11008 6

the file - 625

For manicure work as well as straightforward nail design;
also for delicate work during pedicures



- Speed infinitely adjustable up to 25,000 rpm
- Lightweight, ergonomic, powerful, vibration-free stainless steel handpiece (160 g)
- Professional rotary collet
- 5 LEDs for speed display
- Safety stop

Art. 200620 | EAN 40 43641 20032 9

the file - 620

The entry-level appliance for
sophisticated nail design



Ideal for fingernails and toenails,
whether natural, gel or acrylic.
Minor calluses and hardened skin
can also be removed effortlessly.



- Space-saving, super lightweight casing
- Handpiece in mauve glitter look (75 g)
- Automatic quick-release collet for standard bits (shaft ø 2.332 - 2.350 mm)
- Rest for the handpiece on the appliance
- All-in-one controller (20,000 rpm (clockwise/anti-clockwise rotation, on/off)



- All-in-one controller (clockwise and anticlockwise, on/off)
- Rotary knob with brightly sparkling gemstone (diamond shape)
- Space-saving, elegant design

TECHNICAL DATA
Speed: Infinitely adjustable up to 25,000 rpm · Power supply: 100 - 240 AC, 50/60 Hz, 1,200 mA; Sec. 24 VDC, 1.5 A · Energy consumption at zero load: < 0.3 W · Average efficiency during operation: Level VI



- Rotary knob with brightly sparkling jewel (diamond shape)
- Safety stop (stop in cases of overload)
- Speed display
- Lightest weight appliance on the market

TECHNICAL DATA
Power supply: 100 - 240 VAC, 50/60 Hz, 800 mA; Sec. 24 VDC, 1.5 A · Rotation: 20,000 rpm · Dimensions approx. 110 x 90 x 80 mm (L x W x H) · Weight: Handpiece 75 g

Art. 200630 | EAN 40 43641 20033 6

the file - 620 deluxe
(Manicure/pedicure)

- Including 43-part accessory kit
- Lightest weight appliance on the market!

"The professional class for home use"



- For nail cosmetics and design through to the treatment of callused skin
- Infinitely adjustable up to 20,000 rpm (clockwise/anticlockwise rotation)
- Handpiece in mauve glitter look (only 75 g)
- Separate power supply unit for maximum torque
- All-in-one controller
- Overload protection control
- Rotary knob with brightly sparkling jewel (diamond shape)



Art. 114010 | EAN 40 43641 11006 2

Pedimed (callus remover)
Professional pedicure set for silky soft feet



- Efficiently and gently remove excess skin without any effort
- Quickly and easily removes calluses on your feet and heel from the first application
- Soft and efficient treatment of your skin and nails
- Safe and easy operation due to its ergonomic design
- Slim, lightweight and easy-to-hold file
- Powerful 20,000 rpm (switchable clockwise/anticlockwise rotation)
- High-quality, durable grinding bits with a standard shaft dental diameter
- Also ideal for when travelling
- Including 30 abrasive belts, 5 bits in a designer box, designer case



Art. 143220 | EAN 40 43641 14009 0

Ultra pro S
(manicure/pedicure)

- Including 42-part accessory kit

"Elegant, first class"



- Impressive ease of use
- The speed is infinitely adjustable up to 13,000 rpm using the keys
- The professional clockwise/anticlockwise rotation simplifies work in difficult areas
- An all-rounder for nail care, callus removal and also for nail design
- Practical storage case



Art. 133110 | EAN 40 43641 13019 0

Pedisenso Duo
(manicure/pedicure)

- Including 44-part accessory kit

"The file for diabetics"



- Injury-free working due to "diabetic-friendly safety collet"
- Automatic stop in the case of excessive pressure
- Tried-and-tested powerful technology with planetary gearbox (for callus removal)
- Infinitely adjustable speed control up to 4,000 rpm
- Accessories: See p. 88 for collets
- Practical storage case



Art. 173010 | EAN 40 43641 17001 1

Sensitive
(manicure/pedicure)

Including 42-part accessory kit

“The top class unit
with LED lighting”



- Cold light lamp (LED lighting), safety stop
- Infinitely adjustable speed control up to 12,000 rpm
- Easy to use, dust-protected soft-touch buttons to control the clockwise/anticlockwise rotation
- Power transformer for high-torque power delivery; intensity display with 5 LEDs
- Including practical travel case and microfibre dust cloth

Art. 100030 | EAN 40 43641 10005 6

Emotion
(manicure/pedicure)

Including 41-part accessory kit

“All-round appliance with LED light”



With cold LED spotlight



- For easy and safe care of nails and problem areas such as calloused skin on the ball of the foot or heel
- Speed infinitely adjustable up to 12,000 rpm
- Safety stop prevents possible injuries in the case of excessive pressure
- Convenient button operation for speed and with clockwise/anticlockwise rotation

Art. 155020 | EAN 40 43641 150099

Feeling NEXT GENERATION
(manicure/pedicure)

New: now with USB charging function!

Including 6-part accessory kit

“Battery-operated”



- Shortens and shapes nails, repairs torn (artificial) nails, easily removes cuticles and thickened areas of skin
- New collet for faster bit changes, optimised concentricity
- 2 levels (7,000 rpm, 11,000 rpm)
- Maximum freedom of movement without cables
- Strong and long-lasting battery power



Art. 330212 | EAN 40 43641 33032 3 | White
Art. 230213 | EAN 40 43641 33033 0 | Black

UC-60

Ultrasonic cleaning device 650 ml

NEW!



promed

- Effective cleaning with 45 kHz ultrasonic – thorough and gentle cleaning of jewellery, watches, dentures, nail bits and much more.
- Practical timer with 5 runtimes – selectable programs (90-600 seconds) for flexible application
- 650 ml tank volume – ideal for objects of different sizes

- Compact & lightweight – only about 625 g, space-saving design for home or studio
- Safe & versatile – suitable for water and special cleaning solutions, automatic shutdown after 5 minutes

TECHNICAL DATA

- Dimensions: approx. 198 x 118 x 107 mm • Tank dimensions: approx. 175 x 88 x 49 mm
- Tank capacity: 650 ml • Internal dimensions of cleaning basket: approx. 153 x 52 mm
- Weight: approx. 625 g • Power supply: 220-240 V ~ 50 Hz • Power: 50 W • Ultrasonic frequency: 45 kHz

Art. 330210 | EAN 40 43641 33004 0

UC-50

Ultrasonic cleaning device 550 ml

Used for: Genuine jewellery, glasses, dentures, instruments, grinding bits, CDs/DVDs and much more.

Suitable for tap water



Product video:



Picture only serves as an example.

- Cleans instruments, grinding bits, genuine jewellery, watches, CDs/DVDs, dentures, etc. with ultrasonic, high-frequency vibrations.
- Very stable and quiet
- Automatic switch-off of the UC-50 when the timer expires
- 550 ml tank volume

- Illuminated LCD display; timer function (5 timer settings)
- Also suitable for tap water
- Suitable for disinfectants (improved cleaning and disinfection)
- Accessories: Basket, bit holder, CD holder, watch holder

TECHNICAL DATA

- Dimensions (W x H x D): 200 x 170 x 150 mm • Tank capacity: 550 ml • Weight: approx. 825 g • Power supply: 230 VAC ~ 50/60 Hz • Power: 50 W • Ultrasonic frequency: 43 kHz

promed

promed

Art. 230110 | EAN 40 43641 23001 2 | White
Art. 230115 | EAN 40 43641 23005 0 | Pink

Nailfan

Appliance for the removal of harmful fine particulate

Tested dust protection class "L" – DIN EN ISO 60335-2-69



Product video:



- The Promed nailfan stands out due to its tested dust class "L", its quality, its ease of use and its exceptional performance
- Ergonomic design for fatigue-free working – for you and your customers
- Exchangeable, antistatic upholstered cover
- No cold fingers despite powerful suction
- Low noise with low power consumption
- Very large dust collection area of up to 20 cm



No cold fingers due to the lateral removal of cold air

- Large (11 x 19 cm), easily replaceable fine particulate filter (washable!)
- No filter bag necessary
- Separate dust emptying compartment; antistatic casing
- Including second filter, 1 x Promed dust cloth
- Accessories (p. 91)

TECHNICAL DATA

- Dimensions: (W x H x D) 308 x 140 x 340 mm • Weight: approx. 1.8 kg • Mains power: 230 VAC, 50 Hz • Fan: Power: 22 W, noise: 46 dBA • Max. air volume 180 m³/h

Art. 230120 | EAN 40 43641 23006 7

Nailfan mini

Dust suction

Suitable for fingernails and toenails

Ideal for use in salons and at home



- Exceptionally lightweight appliance that is easy to operate
- 2 adjustable suction levels (low noise)
- High suction capacity on the entire work surface
- Easy replacement and cleaning of the filter (reusable)
- Ergonomic design with customer-friendly working height for the hands



- Exceptionally lightweight appliance that is easy to operate
- 2 adjustable suction levels (low noise)
- Small, handy and transportable – also for travel
- Accessories (p. 91)

TECHNICAL DATA

- Dimensions (W x H x D): 240 x 185 x 60 mm • Weight: approx. 570 g • Mains power: 100 - 240 VAC, 50/60 Hz 0.6 A • Fan: Level I: 3,000 rpm / 24 W, Level II: 4,500 rpm / 36 W

promed

Art. 230121 | EAN 40 43641 230074

Nailfan Mini-X

Professional suction device with battery –
wireless, powerful & mobile

NEW!




- Wireless battery operation – freedom of movement, for use anywhere
- Extremely powerful extraction performance – for thorough dust absorption in everyday studio use
- Maintenance-free brushless suction motor – durable, quiet & dependable

- Up to 8 hours runtime – only 2 hours charging time
- 10 suction stages and 1 high-performance suction mode
- Compact & lightweight – only about 995 g, space-saving design for home or studio
- Accessories (p. 91)

TECHNICAL DATA

• Dimensions: approx. 240 x 225 x 55 mm • Filter dimensions: approx. 194 x 194 mm
• Battery capacity: 4000 mAh • Battery life: 6-8 h • Power supply: 220-240 V ~ 50 Hz • Power: 35 W
• Speed: 4500 rpm • Suction speed: 4.33 m/s • Weight: approx. 995 g



Art. 223000 | EAN 40 43641 220068

UVL 60-X

Professional LED nail lamp with battery –
wireless, powerful & mobile

NEW!



- 42 extra-large high-power LEDs – for maximum light intensity & even results
- Gold-coated LED contacts – improved conductivity & long service life
- Mirror reflectors behind the LEDs to enhance the radiation & optimize curing
- 5200 mAh Battery – Wireless, powerful & long operation without interruption
- Overvoltage protection – max. safety for device & user
- Four timer modes (10/30/60/90 s) – for every application

- High power & low heat function – gentle or extra strong, depending on requirements
- Removable magnetic base plate – ideal for pedicure and quick cleaning
- Operation from top and rear
- Can be mounted under the table
- Robust aluminium handle
- Accessories (p. 91)


TECHNICAL DATA

• Dimensions: approx. 250 x 200 x 90 mm • Number of LEDs: 42 • Weight: approx. 1240 g • Power: 36 W
• Power supply: 220-240 V ~ 50 Hz/DC 12.6 V 2.8 A • Battery capacity: 5200 mAh/57.72 Wh • Wavelength of the UV LED: 365 nm-405 nm • LED life: approx. 50,000 hrs • Timer function: 10 s, 30 s, 60 s, 90 s

Art. 223005 | EAN 40 43641 22004 4

UVL-54 "All-In"

UV LED light curing appliance



Suitable for curing of all UV/LED gel and acrylic materials



- Suitable for fingernails and toenails, for curing of UV + LED gel and acrylic
- Dual power system with 36 LEDs for even curing
- 4 different timer settings (10 s / 30 s / 60 s / 99 s)
- Anti-heat function with 99 s timer
- Removable base (magnetic)



- Automatic sensor system (on/off – integrated light barrier up to 120 s)
- High-definition LEDs: Lifespan approx. 50,000 hours
- Digital display

TECHNICAL DATA
Dimensions (LxWxH) in mm: 221 x 205 x 92; Power supply: 100 - 240 VAC; 50/60 Hz, 0.3 A;
Wavelength of UV LED(λ): 365 nm - 405 nm; LED lifespan: approx. 50,000 h; Tunnel dimensions in mm:
130 x 48 x 127 (W x H x D)

promed

Art. 330040 (white) | EAN 40 43641 33019 4

Art. 330045 (pink) | EAN 40 43641 33020 0

UVL-36 S

UV light for curing of gel and acrylic

Including 4x UV lights
(easily interchangeable)



- Professional nail design for the home
- Ideal for curing gel and acrylic
- 360 degree reflective surface for even curing
- Easy operation due to 4 different timer settings (90 s / 120 s / 30 min. / 60 min.)
- Removable base for easy replacement of the light



- 4 slip-proof rubber feet for greater stability
- See p. 91 for accessories

TECHNICAL DATA
Dimensions: approx. 253 x 221 x 113 mm (LxWxH) • Weight: approx. 1.04 kg • Mains power: 230 VAC,
50 Hz, 36 W • Operating times: 90 s / 120 s / 30 min. / 60 min. • Light socket G23 • Wavelength of the lights:
λ=365 nm

Art. 224010 | EAN 40 43 641 22003 7

LTM-30


LED table magnifying lamp
for comfortable work



- Compact LED table magnifying lamp for your workplace, your workshop, for reading, and much more
- 30 individual LEDs ensure excellent lighting
- 175% magnification while working (3 dpt)
- High-quality 4-inch lens with protective cover
- Easy adjustment of the lens position and height
- Swivelling head and arm joints (multidimensional)



- Ergonomic design and stable base
- Extremely long-life LEDs (high-quality long-life)
- LEDs are much more robust than conventional neon tubes
- Cool shadowless lighting for bright, comfortable working
- Lens size: 4-inch lens (Ø 10.2 cm)
- Colour temperature: 6,200 K (neutral white)
- Dimensions (ØxWxH): 16 x 18 x 39 mm



Art. 330110 | EAN 40 43641 33017 0


LTL-749

LED table lamp
for maximum flexibility



- Elegant and stylish design
- Very bright and even illumination
- No tired eyes, relaxed working
- Ideal for your nail studio, your workplace, on every table, in the bedroom or at school
- 3-level brightness dimming at the touch of a button
- Multidimensional 3-axis adjustment
- 360° all-round adjustable LED tube


- High-quality, long-life LED tube: approx. 40,000 hours
- Energy saving compared with conventional lamps
- LED lamp with integrated glare protection
- Frame made of high-quality aluminium
- Anti-slip base
- Dimensions of the LED tube: approx. 33 cm in length
- Brightness: approx. 5,000 K (corresponds to daylight)





Art. 334000 | EAN 40 43641 33026 2
AC-4000
Air purifier with Wi-Fi function
CADR (clean air delivery rate)
cleaning performance of up to 420 m³/h

Effective against viruses and bacteria,
all kinds of dust, pollen, mites,
smoke/odours



- High filter capacity due to combination filter (HEPA filter H13, activated carbon, photo-catalytic, cold catalytic)
- 3-level aluminium pre-filter (washable)
- UV disinfection by means of ultraviolet light
- 3-stage colour display AQI for air quality (air quality indicator), integrated particle sensor
- CADR (clean air delivery rate) cleaning performance of up to 420 m³/h



- Timer function from 1 to 9 hours
- Temperature indicator and humidity indicator
- The fan has 3 power levels, as well as sleep mode, auto mode, manual mode
- Suitable for rooms up to 70 m²
- Automatic reminder to replace the filter
- Easy to use due to large touch display, child safety lock



Promed AC-4000 – the ideal solution
for clean air in closed rooms.

This high-quality air purifier with combination filter effectively filters viruses, bacteria, as well as all kinds of harmful dust (including fine particulate), mould, pollen, mites, smoke and odours. Up to 99.95% of frequently occurring micro-organisms are removed from the air.

- ✓ Filters particles effectively with an efficiency of up to 99.95%
- ✓ Improved air purification by means of an ion generator: generates 5 million anions/cm³;
- ✓ Improves air quality and promotes good health, particles are safely removed

TECHNICAL DATA

Application area:	approx. up to 70 m ²	Power supply:	220 - 240 VAC, 50 Hz Rated output: 85 W
CADR cleaning performance:	380 - 420 m ³ /h (Clean air delivery rate)	Noise level during operation:	35 dB (A) - 55 dB (A)
Anions concentration:	generates 5 million/cm ³ (improved air quality)	Dimensions in cm:	approx. 41.5 x 25.5 x 69.5
UV lamp:	220 V, 4 W, 254 nm, UV-C, socket T5/G5	Weight:	approx. 8 kg

SCOPE OF DELIVERY

- 1x air purifier AC-4000
- 1x aluminium pre-filter and combination filter
- 1x UVC disinfection tubular lamp
- 1x remote control
- 2x 1.5 V batteries AAA
- 1x instructions for use





Art. 470001 | EAN 40 43641 40085 9 | **white**
Art. 470002 | EAN 40 43641 40086 6 | **black**



Oral irrigator MD-310
Gentle and thorough cleaning for teeth,
gums and interdental spaces

NEW!



- 5 cleaning modes – normal, gentle, pulse, DIY, spot jet
- 310 ml water tank – sufficient for one complete application
- High frequency (1400 pulses/min) – effectively removes plaque and food debris
- 360° rotating nozzle – gets into hard-to-reach places
- Long battery life – up to 35 days on a single charge
- USB-C charging connector – fast and convenient
- IPX7 waterproof – safe to use even in the shower
- Easy & portable – weighs just 270 g, great for on the go
- Accessories/spare parts p. 95

TECHNICAL DATA
• Dimensions: approx. 83.4 × 80.1 × 226.6 mm • Weight: approx. 270 g • Power supply: USB-C charging cable with 5 V/500 mA output • Integrated battery with 3.7 V and 1500 mAh • Tank capacity 310 ml

Health and blood pressure

Blood pressure monitors support better preventative health care and optimise your health.

The first step to better health: monitoring

More and more people are placing a greater focus on monitoring their health – on the one hand, this is due to the better health education among the general population, and on the other, because of the progress in the field of medical research.

Regularly measure and monitor your blood pressure and heart rate with precision. Keep track of your health and identify irregularities in good time.



Art. 402010 | EAN 40 43641 40011 8

PBM-3.5

Upper arm blood pressure monitor



Including cuff, universal (22 - 42 cm), carry case and batteries



WHO Classification
WHO-World Health Organization

- Precise blood pressure measurement on your upper arm
- 90 memory spaces each for 2 users
- Display screen: Systolic, diastolic, heart rate, date, time
- Arrhythmia indicator
- Display of the average value of the last 3 measurements
- Large digital LCD display: approx. 72 x 53 mm

See p. 94 for replacement cuff

- Automatic switch-off function after 3 minutes
- Operates with 4 AA batteries or 6 volt mains adapter
- Includes: Batteries and storage bag
- Measuring range: – blood pressure: 0 - 299 mmHg
– Heart rate: 40 - 199 bpm

Art. 412020 | EAN 40 43641 40059 0

HGP-30

Wrist blood pressure monitor



WHO Classification
WHO-World Health Organization

- 120 memory spaces each for 2 users
- Display: Systolic, diastolic, heart rate, date, time
- Arrhythmia indicator
- Battery level indicator
- Calculation of the average of the last 3 saved values

- Classification of measured values according to the WHO assessment system
- Includes: 2x AAA (1.5 V) batteries, storage box
- Cuff for wrist sizes of 135 - 210 mm
- Measuring range: – Pressure: 30 - 280 mmHg
– Heart rate: 40 - 200 bpm



Art. 421010 | EAN 40 43641 40041 5
INH-2.1R
Ultrasonic inhaler
Gentle treatment for
respiratory diseases



Suitable for allergy sufferers: 100% latex-free
Handy starter set for travel



- High efficiency due to innovative micro membrane nebulisation (mesh technology) with tiny droplets
- Short inhalation time due to very high nebulisation performance (0.2 - 1.0 ml/min.)
- Particularly efficient due to the high respirable aerosol proportion
- Inhalation with mask or mouthpiece – also suitable for children
- Low noise nebulisation
- Small and handy (portable) – ideal for travel
- Filling volume content: max. 8 ml
- Nebulisation performance: approx. 0.2 - 1.0 ml/min.
- Particle size: 5.5 µm
- Size: approx. 42 x 45 x 93 mm
- Includes: Storage bag, mouthpiece, mask for children and adults, connecting piece, batteries (for accessories, see page 94)

ACCESSORIES



Art. 432020 | EAN 40 43641 40078 1

IRT-80

Infra-red thermometer



promed

Precise measurement within seconds



- Contactless infra-red thermometer
- Hygienic and safe measuring
- Measurement of body and object temperature within seconds
- High temperature alarm (visual and acoustic)
- 35 memory spaces
- 2 different display colours
- For normal and high temperature
- Display in °C and °F
- Automatic switch-off (after 10 sec.)
- Measuring range for the body: 32.0 - 42.9 °C
- Measuring range for surfaces/objects: 0 - 100 °C
- Includes storage bag and 2 batteries



Normal temperature display



"High temperature alarm" via colour indication

Art. 431005 | EAN 40 43641 40076 7

PFT-3.7

Digital thermometer

promed

Includes:
Battery and storage box

- Precise body temperature measurement: orally, axillary, rectally
- Acoustic signal when ready to measure
- Acoustic high temperature alarm
- Automatic saving of the last measurement value
- Waterproof
- Automatic switch-off after 10 minutes
- LCD display with digital indicator
- Easy to change batteries
- Measurement range: 32.0 - 42.9 °C

Art. 482020 | EAN 40 43641 40084 2

PM-300 Pro

fingertip pulse oximeter

For measuring the blood oxygen level (SpO₂) and heart rate

NEW!



- OLED display – high contrast, easy to read, even in low light
- SpO₂ & heart rate measurement – fast, accurate results
- Pulse wave indicator – visual control of heartbeat quality
- Compact & lightweight – only about 35 g, ideal for on the go

- Automatic shutdown – energy saving after 8 seconds without signal
- Simple one-button operation – ready for use

TECHNICAL DATA

- Device type: Fingertip pulse oximeter PM-300 Pro • Display type: OLED • Power supply: 1.5 V, 2 x AAA batteries
- Operating voltage: 2.6-3 VDC • Dimensions (mm): 33 x 36 x 58 • Weight: approx. 35 g (without batteries)
- SpO₂ measurement range: 35–100%; Accuracy: 70–100% ±2% (below 70% no requirement); deviation: ±1%
- Heart rate measuring range: up to 250 bpm; Accuracy: ±3 bpm

Art. 481010 | EAN 40 43641 40082 8

PM-200 Pro

Fingertip pulse oximeter

promed

For measuring the blood oxygen levels (SpO₂%) and heart rate (pulse frequency)

Display in portrait and landscape format

including carrying strap



- Display in column or wave form
- Choice of 6 different display formats
- 4 adjustable brightness levels
- Two-colour OLED display
- Battery level indicator

- Automatic switch-off after 8 seconds
- Includes carrying strap and 2 x 1.5 V AA batteries
- Size: approx. 59.4 x 33 x 31.2 mm
- SpO₂: Measurement range: 70% – 100%
- PR: Measurement range: 30 bpm – 254 bpm

Art. M79456 | EAN 40 15588 79456 8


PM 100 connect  Bluetooth® SMART

Pulse oximeter

medisana.



- Data transfer with Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*
- Pulse oximeter for measuring the blood oxygen level (SpO₂%) and heart rate (pulse)
- Two-colour OLED display
- Adjustable display brightness with 10 levels
- Easy "one touch" operation – also ideal for travel



- Oxygen level (SpO₂%) and pulse rate can be displayed in bar or wave form
- 6 different display formats
- Automatic switch-off after 8 seconds
- Battery level indicator
- Includes carrying strap
- Includes batteries (2 x 1.5 V AAA)

Art. M79493 | EAN 40 15588 794933

PM 180

Pulse oximeter

medisana.



- Oxygen level (SpO₂%) and pulse rate can be displayed in bar or wave form
- 4-way rotatable screen display for left and right index fingers
- Automatic switch-off after 8 seconds
- Battery change indicator



- Certified medical device
- Battery level indicator
- Measured values shown on illuminated OLED display
- Single battery operation for sleek design
- Includes: Carrying strap, pouch and battery 1 x 1.5 V (AAA)

Art. M79048 | EAN 40 15588 79048 5

Meditouch 2 connect DUAL

Blood glucose monitor

USB 

medisana.

Including "starter set":
10 test strips, lancing device, 10 lancets,
AST cap, control solution,
case, quick guide, batteries

Only available to order in Germany!



- Precise measurement accuracy according to the new ISO 15197: 2013 due to new test strip technology (GDA-FAD)
- Insufficient fill detection
- Acoustic test reminder
- Memory for 480 measurement results



Starter set

- Measuring unit individually adjustable: mg/dL and mmol/L
- Bluetooth® Smart (4.0) data transfer to the VitaDock+ app for iOS (version 7.0 or higher) and Android (version 4.3 or higher)
- Starter set including lancing device, AST cap, 10 lancets, 10 test strips, control solution, case
- For accessories, e.g. test strips, see p. 95



Art. M79030 | EAN 40 15588 79030 0 | mg/dL

Art. M79034 | EAN 40 15588 79034 8 | mmol/L

Meditouch 2 

Blood glucose monitor

Only approx. 5 seconds measuring time

Including "starter set":
10 test strips, lancing device,
10 lancets, AST cap, control solution,
5 alcohol pads, case, quick guide,
diabetic diary, batteries



- Precise measurement accuracy according to the new ISO 15197: 2013 due to new test strip technology (GDH-FAD)
- Insufficient fill detection
- USB connection for data transfer to VitaDock Online and the VitaDock app for iOS (version 7.0 or higher) and Android (version 4.3 or higher)
- Marking function for measured values before and after meals



Starter set

- Acoustic test reminder
- Average value display for 7/14/30/90 days
- Memory for 480 measurement results
- Automatic ketone warning
- Starter set including lancing device, AST cap, 10 lancets, 10 test strips, control solution, case
- 0.6 µl blood volume
- Large, easily readable display
- For accessories, e.g. test strips, see p. 95



Art. 352430 | EAN 40 43641 35034 5

IT-6 pelvic floor trainer

Therapeutic appliance for training and strengthening the pelvic floor



Including 4 electrodes, 1x vaginal probe, batteries and practical storage case



- Used to prevent or reduce bladder weakness
- Treatment of incontinence after giving birth
- Training the sphincter muscles – for treating stress/urge/mixed and faecal incontinence
- The electrical stimulation helps you regain a feel for the pelvic floor and thus learn to effectively contract the pelvic muscles again



- The appliance fulfils all your needs through the simultaneous treatment of two different regions using 4 electrode connections; the dual outputs mean that two probes can be used at the same time and controlled separately
- For anal and vaginal stimulation, can improve sex life
- See p. 92 for available accessories

Art. 352420 | EAN 40 43641 35030 7

EMT-6 combination appliance

Digital, electrical pain therapy and muscle stimulation appliance

Can also be used in sports



- Appliance for effective pain and muscle treatment
- Designed for the treatment of chronic and acute pain conditions, mobility disorders and for the rehabilitation, maintenance, development and training of muscles and muscle groups (e.g. for sport)
- Can be prescribed by a doctor (reimbursement from health insurance provider)



- Two channel appliance, biphasic asymmetrical, rectangular phase current, direct current neutral
- Scope of delivery: Combined appliance for TENS and EMS in a practical storage case, 4 electrodes (40 x 40 mm), 2 electrode cables, batteries, belt clip
- Various electrodes available, see accessories p. 92

Art. 352410 | EAN 40 43641 35025 3

EMT-4

Digital, electronic pain and muscle stimulation appliance



- 11 pre-set programmes and 2 freely programmable programmes
- Intensity and switch-on protection
- Combination appliance for electronic nerve and muscle stimulation, pain therapy and wellness
- Clock function, standby mode
- Can be prescribed by a doctor (reimbursement from health insurance provider)



- Scope of delivery: Promed EMT-4 appliance, 2 pairs of self-adhesive electrodes (40 x 40 mm), appliance stand, batteries, 2 electrode connection cables, instructions for use, practical storage case, belt pouch
- Various electrodes available, see accessories p. 92

Art. 352020 | EAN 40 43641 35001 7

Tens 1000 s

Analogue, electronic pain therapy appliance



- Separately adjustable frequency and pulse width, choice of 3 amplitude modulations
- Timer: continuously adjustable as well as at 15 and 30 min.
- Can be prescribed by a doctor (reimbursement from health insurance provider)



- Scope of delivery: Promed Tens 1000 s, 4 self-adhesive electrodes (40 x 40 mm), 9-volt battery, 2 contact cables, instructions for use, a practical plastic case
- Various electrodes available, see accessories p. 92





Art. 631010 | EAN 40 43641 60002 0

FB-100

Foot spa



Including 3 pedicure attachments



- Relaxing bubble and vibration massage
- Pleasant red light field
- Keep-warm function
- Massage rollers for foot reflex zone massage
- Massage-supporting insole
- Dry massage function

- 3 exchangeable pedicure attachments: (1) massage attachment, (2) callus remover, (3) brush
- Slip-proof rubber feet
- Comfortable size (suitable up to shoe size 46)
- Power supply: 220 - 240 VAC / 50 - 60 Hz / 80 W

Art. 621010 | EAN 40 43641 60006 8

ITM Pro

hand massage appliance



2 attachments each for finger massage and for acupressure



- Powerful tapping massage
- Infinitely adjustable control of the massage intensity
- Additional red light and warming function
- Various massage attachments: Acupressure and pressure point massage

- Ergonomic shape makes it easy to use
- Non-slip handle
- Extra-long cable
- 230 VAC/ 50 Hz/ 35 W



Art. 340010 | EAN 40 43641 34001 8

"Orchid" douche

Appliance for feminine hygiene



MADE IN GERMANY
10 YEARS WARRANTY

- Ergonomic shape, easy to handle
- Made from stainless steel
- 7 nozzles at the end of the rod
- Including 2 x pressure reducers, 1 x storage box
- Connection thread: G½ inch
- Autoclavable and dishwasher safe
- Made in Germany
- 10-year guarantee

Art. 340020 | EAN 40 43641 34002 5

"Rose" douche

appliance for feminine hygiene



MADE IN GERMANY
2 YEARS WARRANTY

- Ergonomic shape, easy to handle
- Made from non-rusting, high-quality stainless steel
- 10 jets at the end of the rod
- Including pressure regulator and instructions
- Autoclavable and dishwasher safe
- 2-year guarantee
- Connection thread: G½ inch
- Weight: approx. 205 g

Excellent heating products,
excellent quality.



Art. 452120 | EAN 40 43641 400798

KHP-2.3 G

"Sherpa" cosy electric blanket

promed



- 10 temperature levels
- 9 different adjustable times (1-9 hrs.)
- Quick and efficient heating mechanism with 150 watts (Turbo Heat)
- Safety system with overheating protection
- Automatic switch-off according to chosen time setting
- Pleasantly cosy warmth on the entire surface

- Practical controller with illuminated LED control display
- Fluffy, soft upper material, double-sided
- Washable at 40 °C due to removable switch
- Material: 100% polyester flannel
- Size: approx. 180 x 130 cm
- Power: 150 watts

Art. 452110 | EAN 40 43641 40025 5

KHP-2.3

Cosy electric blankets

promed



- 10 temperature levels
- 9 different adjustable times (1-9 hours)
- Quick and efficient heating mechanism with 150 watts (Turbo Heat)
- Safety system with overheating protection
- Automatic switch-off according to chosen time setting
- Pleasant cosy warmth on the entire surface

- Practical controller
- Fluffy, soft upper material
- Washable at 40 °C due to removable switch
- Material: 100% polyester flannel
- Size: approx. 180 x 130 cm
- Power: 150 watts

Art. 452125 | EAN 40 43641 40083 5

KHP-2.3S

Knitted electric blanket

promed



- 6 temperature levels
- Quick and efficient heating mechanism with 160 watts (Turbo Heat)
- Safety system with overheating protection
- Automatic switch-off after approx. 3 hours
- Pleasant cosy warmth on the entire surface
- Practical controller

- Fluffy, soft upper material, double-sided with knit pattern
- Washable at 30 degrees due to removable switch
- Material: 100% polyester
- Size: approx. 180 x 130 cm
- Performance: 160 watts

Art. 612020 | EAN 40 43641 60026 6

HKP-1.8 XXL

Heating pads

promed



- 6 temperature levels
- Heats up quickly with 100 watts (Turbo Heat)
- Automatic switch-off after 90 minutes
- Safety system with overheating protection
- Easy to use, ergonomic operation/switch
- Including elasticated body strap with Velcro fastener and storage bag

- Fluffy, soft upper material
- Washable at 30 °C on a gentle cycle
- Material: 100% polyester
- Size: approx. 60 x 30 cm
- Performance: 100 watts

Art. 614020 | EAN 40 43641 60022 8

NRP-2.0

Shoulder and neck heating pad

promed



- 6 temperature levels
- Automatic switch-off after 90 minutes
- Heats up quickly with 100 watts (Turbo Heat)
- Heats into the neck area
- Safety system with overheating protection
- Fluffy, soft upper material which is comfortable to wear

- Easy to use, ergonomic operation/switch
- Washable at 30 °C
- Material: 100% polyester
- Size: approx. 52 x 56 cm
- Performance: 100 watts

Art. 614010 | EAN 40 43641 600112

NRP-5.4

Back and neck heating pad

promed



- 4 temperature levels
- Automatic switch-off after 90 minutes
- Heats into the neck area
- Safety system with overheating protection
- Washable at 30°C
- Fluffy, soft upper material which is comfortable to wear

- Easy to use, ergonomic operation/switch
- Size when open approx. 60 x 90 cm
- Material: 100% polyester
- Power: 100 watts

Art. 451120 | EAN 40 43641 40028 6

NRP-2.4

Back and neck heating pad

promed



- 10 temperature levels
- 9 different adjustable time settings (10-90 min.)
- Heats up quickly with 100 watts (Turbo Heat)
- Automatic switch-off
- Safety system with overheating protection
- Pleasantly soft removable cover
- Machine washable at 40 °C



- Especially soft, plush fleece cushion
- Perfect hold due to back- and neck-hugging shape
- Practical Velcro strap for individual positioning of the cushion
- Material: 95% plush fleece, 5% nylon
- Size: approx. 39 x 63 cm
- Performance: 100 watts

Art. 450130 | EAN 40 43641 40027 9

HGP-1.7

Heated belt

promed



- 10 temperature levels
- 9 different adjustable time settings (10-90 min.)
- Heats up quickly with 100 watts (Turbo Heat)
- Safety system with overheating protection
- Perfect hold due to adjustable Velcro strap
- Ergonomic operation/switch
- Suitable for an abdominal radius of approx. 125 cm



- Variable wearing comfort: around the stomach or diagonally across the back
- Practical Velcro strap for individual positioning of the cushion
- Washable at 40 °C due to removable switch
- Size: approx. 58 x 30 cm
- Material: 90% plush fleece, 10% nylon
- Performance: 100 watts

Art. 450125 | EAN 40 43641 40062 0

HKP-1.2 H

Cosy heating pad

promed



- 3 temperature levels
- Automatic switch-off after 90 minutes
- Cover and cushion washable at 40 °C due to removable switch
- Safety system with overheating protection
- Easy to use, ergonomic operation/switch
- Heats up quickly with 100 watts (Turbo Heat)



- Including a pad cover (can also be used without cover)
- Fluffy, soft upper material and smooth and soft cover
- Material: 100% polyester
- Size: approx. 40 x 30 cm
- Weight: approx. 0.4 kg
- Power: 100 watts
- Power supply: 230 V/50 Hz

Art. 612015 | EAN 40 43641 60010 5

HKP-1.6 XL

Cosy heating pad

promed



- 6 temperature levels
- Automatic switch-off after 90 minutes
- Heats up quickly with 100 watts (Turbo Heat)
- Easy to use, ergonomic operation/switch
- Includes a 2-sided heading pad cover (can also be used without cover)



- Safety system with overheating protection
- Fluffy, soft upper material and smooth and soft cover
- Cover washable at 30 °C
- Size: approx. 46 x 35 cm
- Material: 100% polyester
- Power: 100 watts

Art. 453010 | EAN 40 43641 40070 5

KFS-100 foot warmer

"No more cold feet"

promed



- 6 temperature levels
- Turbo Heat with 100 watt of power
- Automatic switch-off after 90 minutes
- Safety system with overheating protection
- Ergonomic controller with illuminated LED control display



- Removable, fluffy-soft inner lining
- Inner lining machine-washable at 30 °C
- Exceptionally soft, breathable fabric core
- Large foot space (suitable up to shoe size 46)
- Material: Microfibre (foot warmer), flannel (inner lining)
- Size: approx. 30x30x23 cm

Art. 452210 | EAN 40 43641 40058 3

UBS-2.4 Underbed

Cosy and soothing warmth

promed



- 10 temperature levels
- Rapid heat-up time due to 120 watts (Turbo Heat)
- Safety system with overheating protection
- Automatic switch-off according to chosen time setting (1-9 hrs.)
- Pleasantly gentle and evenly distributed warmth for a restful sleep
- Cosy, breathable material



- Convenient controller with illuminated LED control display
- Machine washable at 40 °C due to removable switch
- Material: 100% cotton (top side) 100% polyester (underside)
- Size: approx. 150 x 70 cm
- Performance: 120 watts



Art. 651025 | EAN 40 43641 60014 3

AL-1300 WS

aroma diffuser with LED with seven different colours



Water tank and nebuliser zone made from antibacterial material – possible bacterial emissions are effectively reduced

With timer function

BPA free and with ultrasonic technology

Dimensions: approx. 15.8 x 22.4 cm



- Aroma diffuser in high-quality wood look and a timeless design, suitable for all room designs; ideal for aromatic oils
- Humidifies and freshens the air in the room for a soothing and therapeutic atmosphere
- Uses ultrasonic technology to create a fine mist of fragrance that spreads throughout the room
- High efficiency due to micro-fine atomisation
- Running time approx. 10 - 14 hours, strong nebulisation
- Suitable for all large and medium-sized rooms
- Automatic switch-off when the tank is empty



- Wellness light in 7 colours, individually selectable
- Low noise and energy-saving
- Timer function: 1/3/6 hrs.
- Includes: Power adapter
- Very large XXL water tank: 1.3 litre
- 100 - 240 VAC / 50 - 60 Hz 24 VDC / 500 mA / 12 W
- See page 54/55 for aromatics



Art. 651030 | EAN 40 43641 60021 1

AL-240 FL

aroma diffuser with flame effect

LED wellness light with flame effect

Dimensions: approx. 210 x 135 x 81 mm



- Aroma diffuser with flame light for aromatic oils
- Humidifies and freshens the air in the room for a soothing and therapeutic atmosphere
- Uses ultrasonic technology to create a fine mist of fragrance that spreads throughout the room
- High efficiency due to micro-fine atomisation
- Automatic switch-off when the tank is empty



- Low noise and energy-saving
- Includes: power adapter and water measuring cup
- Large 240 ml water tank
- 100 - 240 VAC / 50 - 60 Hz 24 VDC / 500 mA / 12 W
- See page 54/55 for aromatics

Art. 651040 | EAN 40 43641 60023 5

AL-400 CL

aroma diffuser with clock and LED with seven different colours

With highlights in noble wood look

Dimensions: approx. 210 x 135 x 81 mm



- Humidifies and freshens the air in the room for a soothing and therapeutic atmosphere – ideal for aromatic oils
- Uses ultrasonic technology to create a fine mist of fragrance that spreads throughout the room
- Wellness light with automatic colour change in 7 colours (colours: can be individually set and dimmed)
- Interval function (30 sec.) of the scented mist
- High efficiency due to micro-fine atomisation



- Automatic switch-off when the tank is empty
- Low noise and energy-saving
- Includes: power adapter and water measuring cup
- Large 400 ml water tank
- 100 - 240 VAC / 50 - 60 Hz 24 VDC / 300 mA / 12 W
- See page 54/55 for aromatics



Art. 651015 | EAN 40 43641 60012 9

AL-150 WS

aroma diffuser with LED in seven different colours

promed

BPA free and with ultrasonic technology

Dimensions: approx. Ø 15.5 x 7.0 cm



- Aroma diffuser in an elegant wood look for aromatic oils
- Humidifies and freshens the air in the room for a soothing and therapeutic atmosphere
- Uses ultrasonic technology to create a fine mist of fragrance that spreads throughout the room
- High efficiency due to micro-fine atomisation
- Automatic switch-off when tank empty
- Wellness light with seven different colours
- Low noise and energy-saving
- Includes: Power adapter
- Economical consumption, large 150 ml water tank
- 100 - 240 VAC / 50 - 60 Hz
24 VDC / 500 mA / 7 W
- See page 54/55 for aromatics



Art. 651010 | EAN 40 43641 60003 7

AL-305

aroma diffuser with LED in five different colours


promed

BPA free and with ultrasonic technology

Dimensions: approx. 15.8 x 22.4 cm



- Aroma diffuser for aromatic oils
- Uses ultrasonic technology to create a fine mist of fragrance that spreads throughout your entire room
- High efficiency due to micro-fine atomisation
- Automatic switch-off when tank empty
- Timer for programming the operating time: 2 and 4 hours
- Wellness light with five different colours
- Low noise and energy-saving
- Includes: Refillable reservoir, power adapter
- 300 ml water tank; runtime: approx. 5 - 9 hrs.
- 100 - 240 VAC / 50 - 60 Hz
24 VDC / 500 mA / 12 W
- See page 54/55 for aromatics



Art. see individual products

Aromatic oil / aromatic essences

(100% essential oil / composition / perfume oil)

promed



Relax aromatic oil (10 ml)

Art. 262010
EAN 40 43641 26005 7



Be Happy aromatic oil (10 ml)

Art. 262011
EAN 40 43641 26006 4



Lavender aromatic essence (10 ml)

Art. 262012
EAN 40 43641 26007 1



Lemon Grass aromatic essence (10 ml)

Art. 262013
EAN 40 43641 26008 8



Orange aromatic essence (10 ml)

Art. 262015
EAN 40 43641 26010 1



Nature-like oil, cinnamon (10 ml)

Art. 262017
EAN 40 43641 26012 5



Eucalyptus aromatic essence (10 ml)

Art. 262018
EAN 40 43641 26013 2

Art. see individual products

Aromatic oil / aromatic essences

(100% essential oil / composition / perfume oil)

promed



Cinnamon-orange aromatic essence (10 ml)

Art. 262019
EAN 4043641 26015 6



Lemon aromatic essence (10 ml)

Art. 262020
EAN 4043641 26016 3



Peach perfume oil (10 ml)

Art. 262021
EAN 4043641 26017 0



Apple-cinnamon perfume oil (10 ml)

Art. 262022
EAN 4043641 26018 7



Wild cherry perfume oil (10 ml)

Art. 262023
EAN 4043641 26019 4



Vanilla perfume oil (10 ml)

Art. 262024
EAN 4043641 26020 0



Rose perfume oil (10 ml)

Art. 262025
EAN 4043641 26021 7



Grapefruit aromatic essence (10 ml)

Art. 262027
EAN 4043641 26023 1



Lime aromatic essence (10 ml)

Art. 262028
EAN 4043641 26024 8

Promed aromatic oils are made in Germany – from the best natural products, processed into 100% essential oil – and provide for a pleasant, sensual fragrance. Treat yourself.

MADE IN GERMANY

- 100% essential oil composition/perfume oil
- Bottle contents: 10 ml
- Packaged in a decorative box
- Made in Germany.

- In 3 languages (DE, EN, FR)
- 26x26x72 mm (LxVxH)
- Suitable for aromatic lamps/aroma diffusers, e.g.: AL-305, AL-150 WS, AL-1300 WS, AL-400 CL, AL-240 FL

Special shoes, soles and accessories

WE MAKE YOU MOBILE!

HOLISTIC QUALITY SHOE CONCEPT

Special shoes for short and long-term care are an essential part of the Promed product range.

The continuous further development in collaboration with orthopaedists, medical specialists, rehabilitation centres, geriatric nurses and diabetes specialists has resulted in one of the most holistic collections for the areas of therapy, rehabilitation and post-operative care.

Medical studies document the high standard of the individual models. The best functionality, high-quality materials and first-class workmanship on the one hand, and innovative ideas and pioneering technology on the other, characterise this holistic, high-quality shoe concept.



Softstep sole concept bandage shoes

In most cases, the outer sole of a bandage shoe doesn't look very spectacular. In most cases, this is not the intention anyway, as the focus is on supporting the healing process. This is precisely why, working together with specialists, Promed has come up with a range of intelligent solutions that can provide both relief for patients and support the healing processes.



Front area without tread

This makes everyday walking easier for people suffering from neuropathy and limping.



Metashank inlay

The special Metashank inlay in the forefoot area relieves the pressure on the metatarsal bones and makes every step easier.




Cylindrical cavities

distributed over the entire sole ensure the best cushioning effect. This provides excellent pressure relief not only in the heel, but also in the mid-foot and forefoot areas. It also helps protect the joints and the spine.




Wider sole

Makes everyday walking easier for people suffering from neuropathy and limping.



Diamond pattern tread

together with the "U-"I" verge, ensures slip resistance indoors and outdoors, for harmonious heel-to-toe motion.



Bevelled heel area

The heel area, which is bevelled by 20°, ensures a safe tread and ideal heel-to-toe motion.



Ergonomically moulded outer sole (foot shaped)

- Suitable for sanding
- Orthopaedically adjustable
- Made of high-quality PU



THE PROMED SOFTSTEP SOLE CONCEPT COMBINES SAFETY AND COMFORT TO GREAT SUCCESS.

Promed – excellent foot temperature

AIRMESH – THE TEMPERATURE-REGULATING UPPER MATERIAL

[Meshline upper material]



AirMesh with EasyDry lining

The 3-layer shaft material (EasyDry, CosyFit, Mesh) ensures that your feet are comfortable and the temperature is regulated, especially in the summer. The EasyDry fibre, which originates from the medical field, transports

moisture away from the foot and dries quickly. Also particularly suitable for allergy sufferers.

Models: TheraLight 1, TheraLight 2-S



AirMesh + Larina new wool

The 3-layer shaft material (new wool, CosyFit, mesh) ensures warm feet without allowing them to overheat, especially in the winter.

Models: Gentlewalk-Hi, Gentlewalk-Lo




The Promed "AirMesh" is particularly breathable due to its mesh structure and the 3 mm open-pore CosyFit lining foam, which makes it temperature-regulating and ensures a great feeling of well-being for your feet.

It originates from the medical field and is free of toxic components. It is particularly suitable for allergy sufferers and diabetics.



COSYFIT – THE BREATHABLE SHAFT MATERIAL

[Masterline insole]



Models:

- all Masterline, apart from München
- all Classicline, apart from München
- Classic 1 and 2.



This material ensures an exceptionally high level of wear comfort and better pressure relief in sensitive areas.

The multi-layer material absorbs the moisture of the foot well and dries quickly after wearing, thereby ensuring a good foot temperature.



In most cases, corns consist of a painful callus with central, deep-pointing (nail-shaped) spores that are caused by chronic pressure on areas where the skin is near the bone.



Pain in the forefoot can have many causes. The most common cause of such pain is inflammation, often due to incorrect footwear and/or overburdening. It can often occur due to foot deformities that result from accidents.



The various types of inflammation often involve joints, tendons, muscle strands or the periosteum. In the case of a plantar fasciitis, the connective tissue in the sole of the foot is affected.



Hammer toe is a toe deformity in which the toe joints contract when they are bent.



A callus is caused by pressure or friction. This results in a slight inflammation of the skin. To protect the underlying cell layers, the skin is locally stimulated so as to thicken the outer epidermis.



Diabetic foot syndrome, also known colloquially as "diabetic foot", is a syndrome which is associated with diabetes, and which occurs most frequently in patients with type 2 diabetes.



Arthritic foot is the term used to describe an inflammation of the joints in the foot or the toes.



The term splayfoot refers to a spreading of the bones of the metatarsus.



Toe deformities are pathological changes in shape that are either congenital or acquired over the course of life. Increased mechanical pressure can lead to painful calluses and inflammation.



Hallux valgus is the medical term for an X-shaped bend to the big toe in the metatarsal joint. The condition almost always occurs together with splayfoot.



The multi-layer **shaft material** absorbs the moisture of the foot well and dries quickly after wearing, thereby ensuring an **excellent foot temperature**.



Particularly suitable for pressure sensitivity of the feet, especially for **diabetes** care and anyone who wants to support healing processes on their feet.



Particularly suitable for rheumatic symptoms such as rheumatic foot changes, claw and hammer toes, hallux rigidus.



The **Softstep** sole concept, with its cylindrical cavities and neuropathic surface (toe area), combines safety and comfort like never before.



The **AirMesh** shaft material is particularly breathable and temperature-regulating thanks to its mesh structure and open-pore CosyFit foam lining. Also particularly suitable for allergy sufferers and diabetics.



The **"washable insole"** principle allows the sole to be changed at any time and washed at 30 degrees. This prevents unpleasant odours.



Natural new **sheep's wool** ensures comfortably warm feet and a wonderfully soft wearing comfort throughout the shoe.



Upper material made from (perforated) **stretch** for especially good breathability. Adapts to almost any foot shape.



Our shoes are manufactured in the EU under strict quality control



This is a mark of the highest quality. The insoles are the substructure of a shoe. The insole creates a durable, dimensionally stable structure, also without an outer sole, which can therefore be provided with different characteristics.



The hard-wearing material, which is visually similar to suede, ensures a healthy breathable temperature in the shoe. It is abrasion-resistant, tear-resistant and easy care.

SYMBOLS



Ideal for orthopaedic use



Washable



Polyethylene velour



SuperComfort sole wider at the forefoot, protection against tripping, slip resistance



Washable at 30 degrees



Cosyfit lining



Flexible foam insert



Micro velour



Terry lining

General information on the range of Promed shoes:



Machine wash at 30°C



Not suitable for tumble dryer



Do not iron



Do not dry-clean



Do not use bleach

The shoes can be disinfected using any commercially available or common household disinfectant.

Overview of shoes and indications



Designation of shoe	Gentle Walk-Hi	Gentle Walk-Lo	München 2 LXL	München 3	Wallgau	Grainau	Flex 2
Picture of the product							
Art. no.	510341	510331	510230	521241 (b.) 521242 (s.)	520210	520220	520230
HiMi no. (Pair)	31.03.03.4058	31.03.03.4059	31.03.03.4059	31.03.03.4151	31.03.03.4152	31.03.03.4149	31.03.03.4150
Page in the 2026 Promed catalogue	67	67	68	68	64	64	63
Main areas of application							
For general problems with the toes and forefoot				✓	✓	✓	✓✓
Ideal for diabetes care	✓✓	✓✓		✓✓	✓✓	✓✓	✓✓
Ideal for rheumatism sufferers	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓
Ideal for neuropathies	✓	✓	✓	✓	✓	✓	✓
Ideal for wheelchair users	✓✓	✓✓	✓	✓✓	✓✓	✓✓	
For almost any foot shape e.g. hallux, hammer toes, arthritis, swollen feet			✓✓		✓	✓	✓
For inflammations, swellings or bandaged feet			✓		✓	✓	✓✓
For foot deformations			✓✓	✓✓	✓✓	✓✓	✓✓
For use after surgery					✓		
For rheumatic feet	✓✓	✓✓		✓	✓	✓	✓
For feet under heavy strain (e.g. obesity)							
For sore feet e.g. traumatic injuries					✓	✓	✓
For pressure sensitivity, inflammations, swellings and lymphatic feet	✓✓	✓✓			✓		✓
Sores in the forefoot area and on the sole of the foot caused by diabetes							
Sores in the forefoot area and on the sole of the foot caused by other indications such as traumatic injuries							
Sores in the rear-foot area and on the sole of the foot caused by diabetes							
Sores in the rear-foot area and on the sole of the foot caused by other indications such as traumatic injuries							

Bandage shoes ✓ = suitable; ✓✓ = particularly suitable



Theramed D1	Theramed D3	Pantomobil D	TheraLight 1	Sanicabrio DXL	Sanicabrio	Sanicabrio LXL	Sanisoft D	Sanisoft
510160	510180	510140	510311	510130	521120	510150	510170	521110
31.03.03.4042	31.03.03.4041	31.03.03.4043	31.03.03.4059	31.03.03.4045	31.03.03.4020	31.03.03.4046	31.03.03.4047	31.03.03.4022
69	69	70	70	65	66	65	66	67
✓✓	✓✓	✓✓	✓✓	✓✓	✓	✓	✓	✓
					✓✓			✓✓
✓✓	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓
✓	✓	✓	✓	✓	✓✓	✓	✓	✓✓
				✓	✓✓	✓	✓	✓✓
✓	✓✓	✓	✓	✓✓		✓✓		
✓✓	✓✓			✓✓		✓✓	✓✓	
✓✓	✓	✓		✓✓		✓✓	✓✓	✓
✓	✓			✓		✓	✓	
					✓✓			✓✓
✓			✓	✓		✓	✓	
✓	✓	✓		✓✓	✓	✓✓	✓	✓✓
✓✓	✓✓	✓✓		✓	✓	✓✓	✓✓	✓
✓✓	✓✓		✓	✓✓		✓✓	✓✓	
				✓✓	✓✓	✓✓	✓✓	✓✓
	✓✓	✓✓		✓		✓	✓	
		✓✓		✓✓			✓	

SHOE DIMENSIONS PROMED MASTERLINE

Model			SaniCabrio DXL, LXL, SaniSoft D			Theramed D 1			Theramed D 3, Pantomobil D			München 2 LXL		
Design			closed toe / closed heel			open toe / closed heel			open toe / open heel			closed shoe		
English sizes *)	German / French sizes	Foot length in mm	Outsole in mm	Insole in mm	Ball size in mm oversize „J“	Outsole in mm	Insole in mm	Ball size in mm oversize „J“	Outsole in mm	Insole in mm	Ball size in mm oversize „J“	Outsole in mm	Insole in mm	Ball size in mm oversize „J“
3½	36	229	253	238	232	249	239	232	243	239	232	243	234	232
4½	37 ⅓	237	256	242	237	253	243	237	249	244	237	253	238	237
5	38	241	264	249	240	260	250	240	253	248	240	256	242	240
6	39 ⅓	250	273	258	245	264	255	245	260	257	245	264	249	245
6½	40	254	277	261	247	269	259	247	264	262	247	273	258	247
7½	41 ⅓	262	285	268	252	277	268	252	273	271	252	277	261	252
8	42	267	295	276	255	285	275	255	280	278	255	285	268	255
9	43 ⅓	275	298	280	260	295	283	260	285	281	260	295	276	260
9½	44	279	307	290	262	298	288	262	295	289	262	298	280	262
10½	45 ⅓	288	312	294	267	307	296	267	298	295	267	307	290	267
11	46	292	315	300	270	312	301	270	307	303	270	312	294	270
12	47 ⅓	301	324	311	275	318	310	275	312	307	275	315	300	275

SHOE DIMENSIONS PROMED MESHLINE

Model			Gentlewalk Hi SaniCabrio Wallgau, Grainau			Gentlewalk Lo SaniSoft München 3			TheraLight I			Flex 2		
Design			closed toe, closed heel			close toe, close heel			open toe, closed heel			closed toe, open heel		
English sizes *)	German / French sizes	Foot length in mm	Outsole in mm	Insole in mm	Ball width in mm comfort width "G"	Outsole in mm	Insole in mm	Ball width in mm comfort width "G"	Outsole in mm	Insole in mm	Ball width in mm comfort width "G"	Outsole in mm	Insole in mm	Ball width in mm Oversize "J"
3½	36	229	247	238	222	247	238	222	247	235	222	247	238	222
4½	37 ⅓	237	254	242	227	254	242	227	250	240	227	254	242	227
5	38	241	263	249	230	259	246	230	259	247	230	263	249	230
6	39 ⅓	250	268	254	235	273	258	235	263	249	235	268	254	235
6½	40	254	279	261	237	279	261	237	273	258	237	279	261	237
7½	41 ⅓	262	284	268	242	284	268	242	281	270	242	284	268	242
8	42	267	293	276	245	291	275	245	291	278	245	293	276	245
9	43 ⅓	275	299	280	250	299	280	250	291	278	250	299	280	250
9½	44	279	306	290	252	306	290	252	291	278	252	306	290	252
10½	45 ⅓	288	311	294	257	311	294	257	311	300	257	311	294	257
11	46	292	320	303	260	313	300	260	313	302	260	320	303	260
12	47 ⅓	301	321	311	265	321	311	265	321	312	265	321	311	265

*) 1 English stitch = ⅓ inch, i.e. 0.84666667 cm



Art. 520230 (in pairs) | HiMi no.: 31.03.03.4145

Flex 2

Unique bandage shoe shaped as slip-on shoe



- Comfortable slip-on shoe for indoor and outdoor use
- Variable heel strap for 3 possibilities to wear
- Improvable fit due to variable heel and instep strap
- Functioning of the heel strap: pivotable or removable with two press studs, heel width adjustable with Velcro fastener
- Broadened forefoot for comfortable fit and secure footing


- Instep strap with Velcro fastener for perfect and infinite adjustment of the instep height
- Ergonomically shaped, orthopaedically adjustable outer sole with neuropathic zone in the forefoot
- Washable at 30°
- Sizes: 36 - 47, comfortable width J
- Colour: black

Art. 520210 (in pairs) | HiMi no.: 31.03.03.4152

Wallgau

Comfortable low bandage shoe
Therapeutic comfort shoe

High-quality,
2-layer insole
for soft
step comfort



Instep Velcro fastener for extra
convenience when putting the
shoe on or taking it off



- Toe box made of comfortable stretch material
- Particularly suitable for pressure-sensitive feet
- Visually attractive upper material made of elastane (stretch) in the forefoot for particularly good breathability; leather look at rear
- Ergonomically shaped, orthopaedically adjustable outer sole with neuropathy zone in forefoot
- Stretch material in the forefoot area allows for a soft and flexible adjustment to almost any foot shape (hallux, hammer toes, arthritis, swollen feet, etc.)
- Ideal for rheumatic feet and for diabetes care
- Washable at 30 °C
- Sizes: 36 - 47, comfortable width G, colour: black




Art. 520220 (in pairs) | HiMi no.: 31.03.03.4149

Grainau

Therapeutic comfort shoe

High-quality,
2-layer insole
for soft
step comfort



Instep Velcro fastener for
extra convenience when
putting the shoe on or
taking it off



- Toe box made of comfortable stretch material
- Particularly suitable for pressure-sensitive feet
- Visually attractive upper material made of elastane (stretch) in the forefoot for particularly good breathability; leather look at rear
- Ergonomically shaped, orthopaedically adjustable outer sole with neuropathy zone in forefoot
- Stretch material in the forefoot area allows for a soft and flexible adjustment to almost any foot shape (hallux, hammer toes, arthritis, swollen feet, etc.)
- Ideal for rheumatic feet and for diabetes care
- Washable at 30 °C
- Sizes: 36 - 47, comfortable width G, colour: black



promed

Art. 510130 (in pairs) | HiMi no.: 31.03.03.4045

Sanicabrio DXL

Over-ankle, extra wide special shoe



- With wide-opening "Cabrio" fastening system at the forefoot and with Velcro fastener at the heel for an easy time putting on the shoe
- Soft and flexible adjustment
- Ideal for swollen feet, hammer toes, for use after operations and with plaster casts and bandaged feet
- Sizes: 36 - 47, extended width J
- Colour: navy blue



Art. 510150 (in pairs) | HiMi no.: 31.03.03.4046

Sanicabrio LXL

Low extra-wide special shoe



- With wide-opening "Cabrio" fastening system at the forefoot.
- Easy to put on thanks to Velcro at the heel
- Low shaft at the ankle
- Ideal for solving foot problems such as bandaged feet, heavily stressed, pressure-sensitive feet, hallux and hammer toes
- Sizes: 36 - 47, extended width J
- Colour: navy blue



Wide-opening "Cabrio" fastening system
at the forefoot and with Velcro fastener at
the heel for an easy time putting on the shoe



promed

Art. 521120 (in pairs) | HiMi no.: 31.03.03.4020

Sanicabrio

Therapy and rehabilitation shoe with Velcro at the forefoot and heel



- Individual adjustment options
- Also suitable as a first aid shoe
- Ideal for extreme pressure sensitivity and severe swelling
- Dipped edge in the sole for stable hold of soft foam insoles, for example
- Terry lining
- Sizes: 36 - 47, comfortable width G
- Colour: navy blue



Wide-opening "Cabrio" fastening system at the forefoot and with Velcro fastener at the heel for an easy time putting on the shoe



Terry lining



promed

Art. 510170 (in pairs) | HiMi no.: 31.03.03.4047

Sanisoft D

Low-level special shoe



- With large volume and various individual adjustment options in the forefoot and instep area
- A temperature-regulating shoe for diabetics and other foot problems
- Easy to put on and take off, and optimum width adjustment
- Especially suitable for wheelchair users
- Sizes: 36 - 47, extended width J
- Colour: navy blue



Easy to slip on due to wide-opening Velcro fastener in the instep and heel area



promed

Art. 521110 (in pairs) | HiMi no.: 31.03.03.4022

Sanisoft

Therapy and rehabilitation shoe in loafer design



- 2 Velcro fasteners in the forefoot area
- Low shaft cut-out under the ankle
- Dipped edge in the sole for stable hold of soft foam insoles, for example
- Terry lining
- Sizes: 36 - 47, comfortable width G
- Colour: navy blue



Velcro fastener



Terry lining

Lightweight and cushioning sole, slip-resistant and slightly wider in the forefoot

Art. 510341 (in pairs) | HiMi no.: 31.03.03.4.058 (Hi)

Art. 510331 (in pairs) | HiMi no.: 31.03.03.4.059 (Lo)

GentleWalk-Hi / GentleWalk-Lo

Fur boots with Airmesh shaft material and new wool lining



- Breathable "AirMesh" upper material with Larina wool lining for a great foot temperature
- 2 wide-opening Velcro fasteners in the forefoot for easy slip-on and a firm hold
- Pressure-free, ideal fit
- Shock-absorbing "SoftStep" outer sole, suitable for asphalt
- Suitable for diabetics and allergy sufferers



(Hi)



(Lo)

- Machine washable at 30 °C
- Sizes: 36 - 47, comfortable width G
- Colour: black

Art. 510230 (in pairs) | HiMi no.: 31.03.03.4054

München 2 LXL

Stretch micro velour loafer for therapy and rehab

promed



- Increased forefoot volume
- 4-fold stretch material at the forefoot allows for soft and flexible adjustment to almost any foot shape (hallux, hammer toes, arthritis, swollen feet, etc.)
- The stretch material follows your every move
- No painful areas of pressure or friction
- Foot temperature friendly model
- Sizes: 36 - 47, extended width J
- Colour: black



The micro-velour shaft offers a soft, individual hold

Inlaid removable insole

Best finish

Increased forefoot volume and 4-fold stretch material in the forefoot area enable a soft and flexible fit for almost any foot shape



Brown – art. 521241 (in pairs) | HiMi no.: 31.03.03.4151
Black – art. 521242 (in pairs) | HiMi no.: 31.03.03.4151

München 3

Low comfort bandage shoe in a stretch-micro velour combination

promed



Centre and heel part made of micro velour for perfect fit

Instep Velcro fastener for extra convenience when putting the shoe on or taking it off

Slightly damping, anti-slip sole

- Instep Velcro fastener for extra convenience when putting the shoe on or taking it off
- Visually attractive upper material made of perforated stretch Dacron for particularly good breathability
- Ergonomically shaped, orthopaedically adjustable outer sole with neuropathic zone in the forefoot
- Particularly suitable for pressure-sensitive feet
- Stretch material in the forefoot area allows for a soft and flexible adjustment to almost any foot shape (hallux, hammer toes, arthritis, swollen feet, etc.)

- Ideal for rheumatic feet and for diabetes care
- Washable at 30 °C
- Sizes: 36 - 47, comfortable width G
- Colour: brown, black

Art. 510160 (in pairs) | HiMi no.: 31.03.03.4042

Theramed D1

Multifunctional 2-Velcro sandal with closed heel

promed




Breathable and easy-care upper material

- With two Velcro fasteners on the forefoot and instep
- Open toe area, instep strap and closed heel area
- Especially recommended for problems in the toe and forefoot area, such as inflammation, swelling and hallux
- Comfortable in the warmer season due to breathable and easy-care upper material
- Sizes: 36 - 47, oversize J
- Colour: navy blue



Can be completely opened with two Velcro fasteners



Art. 510180 (in pairs) | HiMi no.: 31.03.03.4041

Theramed D3

Multifunctional 2-Velcro sandal with open heel

promed



Three Velcro fasteners for individual adjustment in the forefoot, instep and heel areas

- With three Velcro fasteners for individual adjustment in the forefoot, instep and heel areas. Open heel window
- Recommended for pressure sensitivity in the heel area as well as for inflammation and swelling in the heel and toe area
- Temperature-friendly model
- Comfortable in the warmer season due to breathable upper material
- Sizes: 36 - 47, oversize J
- Colour: navy blue



Art. 510140 (in pairs) | HiMi no.: 31.03.03.4043

Pantomobil D

Lightweight, multifunctional, wide-opening special shoe in slipper design



With wide Velcro fastener in the forefoot for a firm hold and individual fit



Breathable cosy fit upper material (three-layer)

- Recommended for pressure sensitivity, swelling or inflammation in the forefoot and/or heel areas
- Foot temperature-friendly model, particularly pleasant in the warmer season due to breathable cosy-fit upper material (three layers)
- Heel strap included
- Sizes: 36 - 47, oversize J
- Colour: navy blue



promed

Art. 510311 (in pairs)

HiMi no.: 31.03.03.4047

TheraLight 1

Rehabilitation and comfort shoe with open toe area and closed heel cap



Putting the shoe on is extremely easy due to Velcro fastener system

Breathable AirMesh upper material with EasyDry lining



- Breathable "AirMesh" upper material with EasyDry lining for an optimal foot temperature
- wide-opening Velcro fastener to make the shoe easy to put on
- Pressure-free, ideal fit
- Shock-absorbing "SoftStep" outer sole, suitable for asphalt
- Suitable for diabetics and allergy sufferers
- Machine washable at 30°
- Sizes: 36 - 47, comfortable width G
- Colour: black

Art. see below (in pairs)

Replacement soft foam insole (double sole)



Replacement insole (double insole) for the Wallgau and Grainau models
 Sizes: 36 - 47, width J
 Art. 550236

- Soft foam insole 1: Cotton EVA (3.6 mm), black
- Soft foam insole 2: Polyester EVA (4 mm), black

Art. see below (in pairs)

Replacement soft foam insole



Sizes: 36 - 47, width J (in pairs)
 Colour: Black
 Art. 550230 (model: Sanicabrio/DXL/LXL/DS, Sanisoft/D)
 Art. 550231 (model: TheraMed D1, TheraLight 1)
 Art. 550232 (model: TheraMed D3, Pantomobil D, TheraLight 2-S)

- High-quality insole with soft step comfort
- Terry-cloth cover
- 4 mm EVA sole

Art. 550235 (in pairs)

Replacement soft foam insole



- High-quality insole with soft step comfort
- Terry-cloth cover
- 4 mm EVA sole
- Suitable for every model in the Promed range
- Sizes: 36 - 47, comfortable width G
- Colour: white

Art. 550240 (in pairs)

Replacement fur insole



suitable for hand washing

- Easy care insoles for sensitive feet
- High-quality EVA insole with 100% sheep's wool covering, 4 mm thick with soft tread comfort
- Suitable for the models: Gentlewalk-Hi, Gentlewalk-Lo
- Sizes: 36 - 47, comfortable width G
- Colour: white



HUBER by **promed**
Fußtechnologie

THE RIGHT FOUNDATIONS MAKE AL THE DIFFERENCE

Art. 552055

Hallux Valgus

To prevent any further bending of the big toe



HUBER
Fußtechnologie




Excellent for the following symptoms:

- Prevention of hallux valgus
- Sharp pain during the heel-to-toe movement phase of the foot
- Dull pain without stress
- Local pressure sensitivity

Particularly suitable for:

- Sufferers of splayfoot and/or a big toe deformity
- Sufferers of hallux valgus



Features:

- Moisture-absorbing, quick-drying, skin-friendly
- Toxin-free soft microfibre material with pleasant wearing comfort

Sizes (1 pair): 35 - 49

Art. 552050

Heel spur

Heel with U-shaped padded structure




Excellent for the following symptoms:

- Sharp pain when putting weight on the foot
- Dull pain without stress
- Local pressure sensitivity

Particularly suitable for:

- Heel spur, stress on the heel spur
- Persons in occupations that require lots of walking or standing on your feet
- Allergy sufferers and diabetics



Features:

- Moisture-absorbing, quick-drying, skin-friendly
- Toxin-free soft microfibre material with pleasant wearing comfort

Sizes (1 pair): **XS**/35-37, **S**/38-40, **M**/41-43, **L**/44-46, **XL**/47-49

Art. 552040 Joint relief sole 1 (up to 80 kg)

Art. 552045 Joint relief sole 2 (over 80 kg)

Joint relief (</> 80 kg)

For targeted, weight-specific pressure relief




Excellent for the following symptoms:


- To relieve stress on the spine and joints
- Complaints regarding the metatarsals (metatarsalgia)
- Burning sole/heel pain, joint problems (ankles, knees, hips)

Particularly suitable for:

- Persons in occupations that require lots of walking or standing on your feet
- Painful, burning, overloaded, tired, pressure-sensitive feet



Helps prevent heavy, tired legs



Helps prevent sweaty feet



Helps prevent burning feet

Features:

- Moisture-absorbing, quick-drying, skin-friendly
- Toxin-free soft microfibre material offering pleasant wearing comfort

Sizes (1 pair):

XS/35-37, **S**/38-40, **M**/41-43, **L**/44-46, **XL**/47-49




Art. 198001 2

Ceramic modelling set

Professional ceramic modelling set for removing, reworking, modelling of UV gel and acrylic, beautifying, polishing. Fast and effective working. Durable material, low heat development

- Contents**
- | | |
|-------------|---|
| Art. 198604 | Ceramic bit, cylinder, rough, green |
| Art. 198603 | Ceramic bit, large cone, fine, red |
| Art. 198606 | Ceramic milling cutter, long cone, medium, blue |
| Art. 198607 | Ceramic milling cutter, pointed, fine, red |
| Art. 198133 | Silicone polisher, pointed, fine
Including rectangular bit box |




Art. 198002 2

Modelling set blue mix

Professional combination set consisting of ceramic and Blue Mix grinding bits. For matting, shortening, removing, reworking and polishing of gel and acrylic. Durable material, low heat development

- Contents**
- | | |
|-------------|---|
| Art. 198480 | Hard alloy bit, cylinder, blue mix series |
| Art. 198481 | Hard alloy bit, large cone, blue mix series |
| Art. 198606 | Ceramic milling tool, long cone, medium, blue |
| Art. 198607 | Ceramic milling tool, pointed, fine, red |
| Art. 198133 | Silicone polisher, pointed, fine
Including rectangular bit box |



Art. 198003 2

Modelling set rainbow

Professional hard alloy bit set with tungsten nitride coating. For removing, modelling, filing, finishing work on artificial nails, polishing; prevents the development of heat.

- Contents**
- | | |
|-------------|---|
| Art. 198491 | Hard alloy bit, large round cone, rainbow series |
| Art. 198490 | Hard alloy bit, cylinder, rainbow series |
| Art. 198492 | Hard alloy bit, narrow cone, rainbow series |
| Art. 198493 | Hard alloy bit, narrow round cone, fine, rainbow series |
| Art. 198133 | Silicone polisher, pointed, fine
Including rectangular bit box |



Art. 198008 2

Modelling set HM-Bitset Pro

Professional high-quality nail modelling set – the “must-have” for every nail designer. Removing, modelling, beautifying and working quickly – now as easy as can be.

- Contents**
- | | |
|-------------|--|
| Art. 198451 | Hard alloy grinder, barrel, combination cut (fine/rough) |
| Art. 198425 | Hard alloy bit, narrow cone, purple series |
| Art. 198442 | Carbide bit, large pointed cone, red series |
| Art. 198557 | Carbide bit, cylinder, cross-dovetailing |
| Art. 198556 | Titanium bit, cylinder, rounded (gold head) |



Art. 198006 2

Master sculpting set

The complete and professional kit for all types of artificial nail modelling.

- Contents**
- | | |
|---------------|---|
| Art. 198410 | Hard alloy grinder, blue series |
| Art. 198421 | Hard alloy grinder, yellow |
| Art. 198450 | Carbide milling tool, blue series |
| Art. 198540 | Diamond milling tool |
| Art. 198031 | Sapphire bit, round |
| Art. 198551/L | Carbide titanium milling tool, (gold head) long |
| Art. 198030 | Sapphire milling cutter, pointed, |
| Art. 961910 | Design hard box, transparent |



Art. 198005 2

Modelling set basic

The perfect basic kit for complete artificial nail modelling.

- Contents**
- | | |
|-------------|--|
| Art. 198421 | Carbide bit, yellow |
| Art. 198450 | Carbide bit, blue series |
| Art. 198540 | Diamond bit |
| Art. 198551 | Carbide titanium bit, (gold head), short |
| Art. 961910 | Design hard box, transparent |



Art. 198014 1

Personal care set

For a complete manicure of natural nails: shortening of nails, cuticle removal, thinning of exceptionally thick nails, elimination of unevenness on the nail plate.

- Contents**
- | | |
|-------------|----------------------------------|
| Art. 190011 | Abrasive disc, sapphire |
| Art. 190067 | Hard skin abrader, sapphire |
| Art. 198540 | Diamond cone, small |
| Art. 198535 | Stone abrader, standard |
| Art. 198541 | Diamond milling cutter, barrel |
| Art. 198031 | Sapphire milling cutter, round |
| Art. 198030 | Sapphire milling cutter, pointed |
| Art. 190041 | Felt cone buffer |



1 3

Art. 198016

Comfort set

For all problem areas on the skin, calluses and nails. Removal of calluses, hardened skin, excess cuticle, nail trimming, matting, smoothing, cleaning and polishing.

Contents	Art. 198535	Stone sander, standard
	Art. 198530	Stone sander, cylinder
	Art. 198564	Abrasive mandrel, plastic
	Art. 198031	Sapphire milling cutter, round
	Art. 198543	Diamond milling cutter, small, pointed
	Art. 198034	Sapphire bit, long
	Art. 198540	Diamond bit
	Art. 198541	Diamond milling cutter, barrel
	Art. 190086	Goat hair brush
	Art. 190041	Felt cone buffer
	Art. 190021	Abrasive cone, sapphire
	Art. 190011	Abrasive disc, sapphire
	Art. 198562	Abrasive band, medium (4 pcs.)
	Art. 961910	Design hard box, transparent

Art. 198009

Natural set

For a complete natural nail manicure: Shortening of nails, cuticle removal, thinning of exceptionally thick nails, elimination of unevenness on the nail plate, through to the perfect finish with shiny polishing.

Contents	Art. 190011	Abrasive disc, sapphire
	Art. 198020	Abrasive cone, diamond
	Art. 198540	Diamond cone, small
	Art. 198535	Stone abrader, standard
	Art. 198130	Silicone polisher, fine, yellow
	Art. 198131	Silicone polisher, medium, black
	Art. 198132	Silicone polisher, rough, green
	Art. 190041	Felt cone buffer
	Art. 198571	Soft polisher, textile

Art. 198010

Callus set

For the removal of very hard calluses and hardened skin from the hands and feet and other parts of the body, such as the elbows; also for cuticle removal.

Contents	Art. 190021	Abrasive cone, sapphire
	Art. 190067	Hard skin abrader, sapphire
	Art. 198540	Diamond cone, small

* Autoclavable (sterilisation of bits at a temperature of between 121 - 134 °C at a pressure of approx. 2 bar)

1 = for treatment of natural nails 2 = for treatment of artificial nails 3 = for skin treatment

Note: All Promed grinding bits have a shaft diameter of a standard dental diameter of 2.332 to 2.350 mm.

Ceramic bits

More lightweight than carbide bits, durable, efficient

Suitable for removing gel and acrylic, filing, engraving, smoothing, polishing and reshaping as well as removing dead skin. Extremely good treatment of nails possible, smooth surfaces, pleasant feel on the nail, fast, easy handling, good running properties due to reduced vibration, very corrosion-resistant, low heat development.

Art. 198600

Ceramic milling cutter, large cone, fine, red *

Suitable for reworking very thickened nails and calluses as well as artificial nails; smooth running properties, extremely hard, durable, sharp. Speed range: up to 30,000 rpm

Art. 198601

Ceramic milling cutter, large cone, medium, blue *

Suitable for treating very hard, thickened nails and calluses as well as artificial nails. Ideal for removing UV gel and acrylic effectively. Also suitable for reworking. Lower weight compared with carbide bits, smooth running properties, extremely hard, durable, sharp. Speed range: up to 30,000 rpm

Art. 198602

Ceramic milling cutter, large cone, rough, green *

Suitable for treating very hard, thickened nails and calluses as well as artificial nails. Ideal for removing UV gel and acrylic effectively. Lower weight compared with carbide bits, smooth running properties, extremely hard, durable, sharp. Speed range: up to 30,000 rpm

Art. 198604

Ceramic bit, cylinder, rough, green *

Matting of the tips, shortening of artificial nails, removal of gel and acrylic. Very good cutting performance, pleasant feel on the nail. Speed range: up to 20,000 rpm

Art. 198603

Ceramic bit, large cone, fine, red *

An absolute "all-rounder" in its field. The unique cone shape allows you to remove gel and acrylic (also 3-phase gels) with little risk. It also provides an optimum basis for further nail modelling work. Very good cutting ability. Speed range: up to 25,000 rpm

Art. 198605

Ceramic milling cutter, long cone, fine, red *

For the finishing, gentle polishing and reshaping of your nails. Speed range: up to 20,000 rpm

Art. 198606

Ceramic milling cutter, long cone, medium, blue *

For removing material, for straightforward surface treatment, possibly also for shortening irregular nails. Speed range: up to 20,000 rpm

Art. 198607

Ceramic milling cutter, pointed, medium, red *

Narrow bit for work at the nail fold/edge, easy trimming or work below the nails. Speed range: up to 35,000 rpm

Art. 198608

Ceramic milling cutter, pointed, medium, blue *

Narrow bit for work at small and narrow spots, within easy reach without scratching the nails. For reworking modelled nails, for removing residues below the nails. Speed range: up to 35,000 rpm

Blue mix bits

Very durable material, prevents the development of heat

Blue mix series: Carbide bits with aluminium nitride coating Particularly rapid and speedy treatment of artificial nails.



Art. 198480

Hard alloy bit (with aluminium nitride coating), cylinder, blue mix series *

Matting of the tips, shortening of artificial nails, removal of gel and acrylic. Coating prevents the development of heat, very durable material. Speed range: up to 25,000 rpm

2



Art. 198481

Hard alloy bit (with aluminium nitride coating), large cone, blue mix series *

An absolute "all-rounder" in its field. The unique cone shape allows you to remove gel and acrylic (also 3-phase gels) free of risk. It also provides an optimum basis for further nail modelling work. Coating prevents the development of heat, very durable material. Speed range: up to 25,000 rpm

2

Rainbow bits

Very durable material, prevents the development of heat

Rainbow series: hard alloy bits with tungsten nitride coating Particularly rapid and speedy treatment of artificial nails.



Art. 198491

Hard alloy bit (with tungsten nitride coating), large round cone *

An absolute "all-rounder" in its field. The unique cone shape allows you to remove gel and acrylic (also 3-phase gels) free of risk. It also provides an optimum basis for further nail modelling work. Tungsten nitride coating prevents the development of heat, very durable material. Speed range: up to 25,000 rpm

2



Art. 198490

Hard alloy bit (with tungsten nitride coating), cylinder *

Matting of the tips, shortening of artificial nails, removal of gel and acrylic. Tungsten nitride coating prevents the development of heat, very durable material. Speed range: up to 25,000 rpm

2



Art. 198492

Hard alloy bit (with tungsten nitride coating), narrow cone, fine *

Gets acrylic and (building) gel into shape. Due to its size, it is especially suitable for finishing work. Tungsten nitride coating prevents the development of heat, very durable material. Speed range: up to 30,000 rpm

2



Art. 198493

Hard alloy bit (with tungsten nitride coating), narrow small round cone, fine *

Shapes acrylic and (build-up) gel and files these down very well. Its shape matches the shape of the nail perfectly and allows for greater work results. Tungsten nitride coating prevents the development of heat, very durable material. Speed range: up to 30,000 rpm

2

* Autoclavable (sterilisation of bits at a temperature of between 121 - 134 °C at a pressure of approx. 2 bar)

1 = for treatment of natural nails 2 = for treatment of artificial nails 3 = for skin treatment

Note: All Promed grinding bits have a shaft diameter of a standard dental diameter of 2.332 to 2.350 mm.

Titanium bits / carbide bits

for extreme loads, long service life

These high-quality bits made of a titanium-carbide alloy do not heat up as quickly, which means they are more pleasant for your customers during use.



Art. 198551

Titanium bit, cylinder, short (gold head) *

Matting of tips, shortening of artificial nails, removal of gel and acrylic material, very durable. Speed range: up to 25,000 rpm

2




Art. 198551/L

Titanium bit, cylinder, long (gold head) *

Matting of tips, shortening of artificial nails, removal of gel and acrylic material, very durable. Speed range: up to 25,000 rpm

2




Art. 198556

Titanium bit, cylinder, rounded (gold head) *

Filing of the tip transition to the natural nail, blending of the tips, shortening of artificial nails, thin filing of the nails' edge; removal of gel and acrylic, speed range: up to 25,000 rpm

2



Art. 198557

Carbide bit, cylinder, cross-dovetailing *

Special bit for both running directions; rapid removal of modelling material (gel and acrylic), special coating prevents heat build-up; very fine polish on the nail, very durable, speed range: up to 25,000 rpm

2

Hard alloy bits

very durable material

Hard alloy bits allow particularly fast and effective work on artificial nails. They must not be used on natural nails. Typical areas of use: Filing of the tip transition, thin filing of the surface, filing of the underside of the nail, shortening, countersinking of rhinestones.

Hard alloy bits, purple series, drop – fine



Art. 198455

Carbide bit, drop *

Nail and skin smoothing as well as nail art work. Speed range: optimum at 12,000 rpm, max. up to 20,000 rpm

2




Art. 198435

Hard alloy bit, long cone *

Carbide bit enables a wide range of applications: Cleaning of the underside of the nail at the nails' edge. Removal of gel or acrylic residues and carrying out of finishing work. Speed range: up to 40,000 rpm

2



Art. 198425

Hard alloy bit, narrow cone *

Particularly effective with "lifting". Finishing work and quick removal of gel or acrylic, also near the nail fold. Speed range: up to 30,000 rpm

2




Art. 198445

Hard alloy bit, large cone *

An absolute "all-rounder" in its field. The unique cone shape allows for the risk-free removal of gel and acrylic (including 3-phase gels) without generating heat. It also provides an optimum basis for further nail modelling work. Extremely durable material. Speed range: up to 20,000 rpm

2



Particularly suitable
for 3-phase gels >

For fine
treatments only >

Hard alloy bits, blue series, cross-dovetailing – rough

Art. 198410 2
Hard alloy grinder, round *
Particularly suitable for countersinking jewels in artificial nails. Speed range: up to 40,000 rpm

Art. 198420 2
Hard alloy grinder, cone *
Shapes acrylic and (build-up) gel and files these down very well. Its shape matches the shape of the nail perfectly and allows for greater work results. No heat build-up.
Speed range: up to 20,000 rpm

Art. 198430 2
Hard alloy grinder, small cone *
Gets acrylic and (building) gel into shape. Its size is particularly suitable for finishing work.
Speed range: up to 30,000 rpm

Art. 198450 2
Hard alloy bit, rounded cylinder *
Very quick and rough removal of the modelling (acrylic and gel). Very low heat generation.
Speed range: up to 20,000 rpm

Hard alloy bits, red series, cross cut – medium

Art. no. 198422 2
Carbide bit, cone *
Optimum removal of all types of gels - particularly suitable for refilling. The cone shape is particularly suitable for finishing the modelling, such as the shaping work.
Speed range: up to 30,000 rpm

Art. 198442 2
Hard alloy grinder, bigfootball *
Due to its shape, the small pointed cone is ideal for the precise removal of gel and acrylic, also in difficult areas just before the nail fold. Speed range: up to 25,000 rpm

Art. 198443 2
Hard alloy grinder, small football *
Due to its shape, the small pointed cone is ideal for the precise removal of gel and acrylic, also in difficult areas just before the nail fold. Speed range: up to 30,000 rpm

Art. 198452 2
Hard alloy grinder, rounded barrel *
The rounded cylinder shape enables the safe and even removal of gel and acrylic over the whole nail plate without generating heat. Speed range: up to 20,000 rpm

Hard alloy bits, yellow series, cross cut – extra fine

Art. 198434 2
Hard alloy grinder, long *
Carbide bit. Its fine surface structure (cross-cut) cleans the underside of the nails perfectly. Removal of acrylic or gel residues, finishing work.
Speed range: up to 40,000 rpm

Art. 198421 2
Carbide bit, cone *
Carbide bit. For reworking the modelling and the removal of residues.
For the perfect cleaning of the underside of the nail. Speed range: up to 30,000 rpm

Art. 198461 2
Hard alloy grinder, flame *
Due to its shape, the flame carbide bit can access any area of the nail, no matter how difficult. For removing adhesive or material underneath acrylic and gel nails. Can also be used to clean the natural nail after regrowth. Speed range: up to 40,000 rpm



Hard alloy bits, special bits

Art. 198451 2
Hard alloy grinder, barrel, combination cut (fine/rough) *
Combination cut (fine + rough elements), very durable material, rapid removal without heat development; very smooth surface structure on the nail.
Speed range: up to 20,000 rpm

Art. 198453 2
Hard-alloy grinder, narrow cone, combination cut (fine/rough), orange series
For finishing the modelling and removing residues, finishing work and quick removal of gel or acrylic without heating, a very smooth surface structure on the nail, very durable material, can also be used for pedicures the removal of calluses. Speed range: up to 20,000 rpm.

Art. 198423 2
Carbide bit, left-handers, cone, cross-dovetailing, medium *
Especially suitable for left-handed operation; quick removal of modelling (gel and acrylic) fine surface texture on the nail; very durable material.
Speed range: up to 20,000 rpm

Art. 198470 2
Special milling cutter with cross-cut, green series, rough
For the treatment of heavily thickened, gryptic nails and calluses.
Speed range: optimum at 12,000 rpm, max. up to 20,000 rpm

Art. 198471 2
Special milling cutter with cross-cut, orange series, fine
Removal of hardened skin and thick calluses as well as nail art work. With external cooling channels, for cool working in the dry and wet technique.
Speed range: optimum at 10,000 rpm, max. up to 15,000 rpm

Steel milling cutter – stainless

Stainless steel, hygienic, durable – for precise work in pedicures/podiatry, ideal for working on corns, calluses, hardened skin and in the nail folds

Art. 197005 3
Round cutter with cross-cut (2.3 mm), medium *
For removing corns and hardened areas of skin. Easy-cut toothing with cross-cut. Speed range: optimum at 15,000 rpm, max. up to 30,000 rpm

Art. 197015 1 3
Fissure cutter (1.2 mm), medium cut *
For work on the nail fold. With cross-cut.
Speed range: optimum at 15,000 rpm, max. up to 50,000 rpm

Art. 197021 3
Hollow cutter (trepan, 2.3 mm) *
For the removal of calluses and hardened skin. Cutting edge with toothing.
Speed range: optimum at 4,000 rpm, max. up to 6,000 rpm

Art. 197025 3
Rotary tissue punch, ground (2.7 mm) *
For the gentle removal of deeper corns/calluses With ground cutting edge.
Speed range: optimum at 4,000 rpm, max. up to 6,000 rpm

Art. 197030 1 2 3
Nail cleaner, short *
For the gentle removal of excess cuticle and nail fold deposits. Hexagonal working part without sharp cutting edges.
Speed range: optimum at 5,000 rpm, max. up to 30,000 rpm

Art. 197031 1 2 3
Nail cleaner, long *
For the gentle removal of excess cuticle and nail fold deposits. Hexagonal working part without sharp cutting edges.
Speed range: optimum at 5,000 rpm, max. up to 30,000 rpm

Diamond bits

Fine diamond particles on the grinding bit allow for hygienic and quick skin and nail treatment. Great durability and longevity with continuous use. Smoothing, cleaning, removal of adhesive residue, also in the nail fold. Safe working in difficult areas.

- Art. 198505

1

2

Nail cutter (only top is diamond-coated) *

Ideal for shortening nails. Adjacent skin is not damaged. Safety bevel on the front edge.

Speed range: optimum at 15,000 rpm, max. up to 20,000 rpm
- Art. 198502

3

Hard skin abrader, cone (super coarse grit) *

Removal of hardened skin and thick calluses. With external cooling channels, for cool working in the dry and wet technique.

Speed range: optimum at 10,000 rpm, max. up to 15,000 rpm
- Art. 198501

3

Hard skin abrader, conical (super coarse grit) *

Removal of hardened skin and thick calluses.

With external cooling channels, for cool working in the dry and wet technique.

Speed range: optimum at 10,000 rpm, max. up to 15,000 rpm
- Art. 198585

3

Diamond milling cutter, rounded barrel, big, rough, hard skin abrader (especially for wet technique) *

For professional pedicure. No heat build-up. For treating uneven skin and calluses.

Suitable for use on dry or damp skin. Speed range: up to 20,000 rpm
- Art. 198586

3

Diamond milling cutter, rounded barrel, small, rough, hard skin abrader (especially for wet technique) *

For professional pedicure. No heat build-up. For treating very thick calluses. Suitable for use on both dry and damp skin.

Speed range: up to 20,000 rpm
- Art. 198020

1

3

Abrasive cone, diamond *

For the sanding of uneven, grooved nail surfaces. Speed range: up to 20,000 rpm
- Art. 198540

1

3

Diamond bit, small cone *

For cuticle removal, the matting of transitions and the natural nail, treatment of difficult edges. Speed range: up to 20,000 rpm
- Art. 198541

1

3

Diamond milling cutter, barrel *

For thinning thickened nails and treating hard corn edges.

Speed range: up to 20,000 rpm
- Art. 198542

1

3

Diamond milling cutter, football *

For the sanding of fungal nails; traces of nail fungus can be easily removed this way.

For cuticle removal. Speed range: up to 20,000 rpm
- Art. 198543

1

3

Diamond milling cutter, small football *

For grinding fungal nails; traces of nail fungus can be easily removed this way.

For cuticle removal. Speed range: up to 20,000 rpm

Diamond bits with cooling channels and carbon shaft

Diamond bits with a high diamond particle density and an integrated cooling channel enable an exceptional removal performance. The cooling channel reduces the heat generation and at the same time ensures longer operating times through the continuous removal of the unwanted material. Suitable for dry and wet technique.

- Art. 198544

1

3

Diamond bit with cooling channels, pointed, small

With cooling channels, ideal for finishing and smoothing of calluses and nails in sensitive areas.
- Art. 198545

1

3

Diamond bit with cooling channels, round

With cooling channels, ideal for the gentle finishing and smoothing of hardened skin and corns.
- Art. 198546

1

3

Diamond bit with cooling channels, thin, rounded

With cooling channels, ideal for smoothing and finishing calluses and nails. Allows for work at lower temperatures.
- Art. 198547

1

3

Diamond bit with cooling channels, conical, thick

With cooling channels, for precise thinning and finishing of hardened skin and nails.

Allows for work at lower temperatures.

Sapphire bits

Sapphire is “softer” than diamond and ideal for the care and treatment of natural nails, thickened nails and calluses or other foot and nail problems.

- Art. 198030

1

3

Sapphire milling cutter, pointed

For treating difficult corners and folds, careful cuticle removal.

Speed range: up to 25,000 rpm
- Art. 198031

1

3

Sapphire bit, round

For the cleaning of the nail fold and the removal of calluses.

For the grinding of ingrown nails. Speed range: up to 25,000 rpm
- Art. 198034

1

3

Sapphire milling cutter, long

For treating the underside of the nail, e.g. removing dirt or traces of nail fungus.

Speed range: up to 25,000 rpm
- Art. 190011

1

Abrasive disc, sapphire

For shortening and shaping fingernails and toenails. Speed range: up to 10,000 rpm
- Art. 190027 fine

1

3

Art. 190029 medium

Art. 190031 rough

Abrasive cone, sapphire, small *

For removing hardened skin or calluses on the sole of the foot. Can also be used to remove corns. Speed range: up to 20,000 rpm

* Autoclavable (sterilisation of bits at a temperature of between 121 - 134 °C at a pressure of approx. 2 bar)
1 = for treatment of natural nails 2 = for treatment of artificial nails 3 = for skin treatment

Note: All Promed grinding bits have a shaft diameter of a standard dental diameter of 2.332 to 2.350 mm.

Sapphire bits

Sapphire is “softer” than diamond and ideal for the care and treatment of natural nails, thickened nails and calluses or other foot and nail problems.

- Art. 190024 fine
- Art. 190026 medium
- Art. 190028 rough
- Art. 190021 rough (not autoclavable)

Abrasive cone, sapphire *

For the grinding of hardened skin, calluses and very woody nails.

Speed range: up to 20,000 rpm

- Art. 190066 fine
- Art. 190064 medium
- Art. 190063 rough
- Art. Art. 190067 rough (not autoclavable)

Hard skin abrader, sapphire *

Particularly suitable for removing large and thick calluses such as on the heel or ball of the foot. Speed range: up to 20,000 rpm

Stone bits

The stone grinding bits are especially suitable for thinning out thicker nails.

- Art. 198530
- Stone abrader, cylinder

For smoothing unsightly nail plates, for thickened, woody nails.

Speed range: up to 12,000 rpm

- Art. 198535
- Stone abrader, standard

For smoothing and matting the nail plate. To prepare for the application of artificial nails.

Speed range: up to 12,000 rpm

Mandrels – abrasive cap holder

For the attachment of abrasive bands for modelling work and callus removal. Easy replacement of the abrasive caps due to the spring action of the cap holder. The new mandrel version effectively compensates for the tolerances of different abrasive caps.

- Art. 198583
- Abrasive mandrel, metal

Suitable for attaching abrasive belts for work on the modelling of artificial fingernails.

Speed range: up to 25,000 rpm

- Art. 198564
- Abrasive mandrel, plastic *

Suitable for attaching abrasive belts for work on the modelling of artificial fingernails.

Speed range: up to 25,000 rpm

- Art. 198565/L
- Abrasive mandrel, plastic, long shaft

(for 4030/S/SX/SX2, 4040, 5030S, 5040-SX, 5040-SX2) *

Suitable for attaching abrasive belts for work on the modelling of artificial fingernails.

Speed range: up to 25,000 rpm

Sandpaper bits

Sanding cylinders are available in various grain sizes. The cylinders are pulled over a mandrel, thereby turning the trusty hand file into a neat, well-round tool.

- Art. 198561 fine

Abrasive band, fine, 100 pcs.

Suitable for sanding the surface of the nail and for delicate work.

- Art. 198562 medium

Abrasive band, medium, 100 pcs.

Particularly suitable for removing gel.

- Art. 198563 rough

Abrasive band, rough, 100 pcs.

Suitable for removing acrylic materials.

- Art. 198590 fine

Abrasive band ZEBRA, fine, 100 pcs.

With “anti-dust” coating. Suitable for sanding the surface of the nail and for delicate work.

- Art. 198591 medium

Abrasive band ZEBRA, medium, 100 pcs.

With “anti-dust” coating. Particularly suitable for removing gel.

- Art. 198592 rough

Abrasive band ZEBRA, rough, 100 pcs.

With “anti-dust” coating. Suitable for removing acrylic materials.

- Art. 198587

Sandpaper cap (10 mm), medium, 100 pcs.

Professional sandpaper cap for pedicures. For the removal of strong and thick calluses and hardened skin on the ball of the foot and heel.

- Art. 198588

Sandpaper cap (13 mm), medium, 100 pcs.

Professional sandpaper cap for foot care. For the removal of strong and thick calluses and hardened skin on the ball of the foot and heel.

- Art. 198581

Abrasive mandrel, rubber (10 mm for sandpaper cap)

Abrasive mandrel to fit sandpaper caps for the treatment of calluses.

Speed range: up to 15,000 rpm

- Art. 198580

Sanding mandrel, rubber (13 mm for sandpaper cap)

Abrasive mandrel to fit sandpaper caps for the treatment of calluses.

Speed range: up to 15,000 rpm



Polishing bits

Polishing, shine removal, shine creation, matting, smoothing unevenness in natural nails and much more. The various polishing bits, from the silicone polisher to the felt bit through to the textile soft polisher, offer a wide range of applications in nail care




Art. 190041
Felt cone buffer
For smooth nails and nail edges to give a shiny finish to the treatment.
Speed range: up to 10,000 rpm

1 2



Art. 198133
Silicone polisher, pointed, fine
For polishing, shining and matting. Speed range: up to 8,000 rpm

1 2




Art. 198130
Silicone polisher, fine, yellow
For polishing, shining and matting.

1 2



Art. 198131
Silicone polisher, medium, black
For polishing, shining and matting.

1 2



Art. 198132
Silicone polisher, rough, green
For polishing, shining and matting. Removal of minor unevenness in natural nails. **Buffer replacement!**

1 2

Note:

Three different grain sizes for the fine polishing of natural and artificial nails.
Speed range: up to 8,000 rpm.
(Warning: If used at speeds over 8,000 rpm, the nails discolour with the colour of the silicone polisher)



Art. 198572
Soft polisher, textile, including polish
Textile including polish creates a uniquely shiny finish on natural nails, as the polishing paste is already integrated into this polishing bit. Speed range: up to 10,000 rpm.

1 2



Art. 198570
Soft polisher, leather
For polishing natural nails. Also for artificial nails when combined with a polishing paste. The surface structure of the leather creates a beautiful lustre on the plate of natural nails. Speed range: up to 10,000 rpm

1 2




Art. 198571
Soft polisher, textile
For polishing artificial and natural nails to a shine in combination with a polishing paste. Speed range: up to 10,000 rpm

1 2




Art. 190086
Polishing brush, goat hair
For cleaning nails and material. Can be used to remove dust residues from nails in hard-to-reach areas, such as around the cuticles. Made from soft goat hair. Speed range: up to 8,000 rpm

2



*Autoclavable (sterilisation of bits at a temperature of between 121 - 134 °C at a pressure of approx. 2 bar)
1 = for treatment of natural nails 2 = for treatment of artificial nails 3 = for skin treatment

Note: All Promed grinding bits have a shaft diameter of a standard dental diameter of 2.332 to 2.350 mm.



Useful professional accessories

Art. 198915
Bit holder, plastic
Space for 13 bits (supplied without bits), white, 100 x 30 mm

Art. 198900
Acrylic bit holder
Space for 10 bits (supplied without bits), transparent, 95 x 50 mm

Art. 198910
Cleaning brush for bits, brass
For cleaning grinding bits.

Art. 198920
Cleaning brush for bits, nylon
For cleaning of drill bits. Particularly gentle on the bit blades.

Art. 198930
Cleaning brush for bits, nylon
For cleaning of drill bits. Particularly gentle on the bit blades.

Art. 940602 (10 pieces) / Art. 940601 (1 piece)
Promed towel
white, small, 30 x 50 cm, scope of delivery: Pack of 10 or single

Art. 952010
Manual for professional nail design
This small, illustrated book discusses the most important techniques for the use of electrical files and covers all the key topics.
68 pages, 80 illustrations, ISBN 3-9809164-0-5, in German

Art. 190051 – white
Gripping collet attachment – white (for Pedisenso Duo)
fixed screw collet

Art. 190052 – yellow
Collet – yellow (for Pedisenso Duo)
with safety slip disc for diabetics

Handpieces

(Also separately available)



Art. 200601
Promed handpiece KR6
separate for Promed 620



Art. 292001
Promed handpiece PS7
separate only for Promed 525 / 625 / 725 / 825
with locking mechanism (not for 520 / 620)



Art. 201004
Promed handpiece MR10
separate for Promed 1030 with locking mechanism



Art. 202003
Promed handpiece MR20
separate for Promed 2020 (from 2015) with locking mechanism



Art. 203031
Promed handpiece The File 3030-X
separate for Promed 3030-X, white stainless steel,
(not available separately)

NEW



Art. 204002
Promed handpiece POM
separate for Promed 4030-SX2, white



Art. 204007
Promed handpiece 4040
separate for Promed 4040, white, with drawer for filter bags
(not available separately)

Accessories, foot pedals

Art. 290011

NEW

Foot pedal 825
ON/OFF foot pedal with robust plastic housing (100 x 70 x 30 mm) and rubber feet. Colour: black, connection cable with 5-pin plug, cable length: 150 cm. Connection to devices **825**

Art. 290004

Foot pedal 1
ON/OFF foot pedal with robust plastic housing (100 x 70 x 30 mm) and rubber feet. Colour: black, connection cable with 5-pin plug, cable length: 150 cm. Connects to appliance models **1030, 2020, 2030** and **3020**

Art. 290007

Foot pedal 5
Same as foot pedal 1, ON/OFF foot pedal with plastic casing. Only connects to the models: **4030-SX, 4030-SX2**

Art. 290009

Foot pedal 6
Same as foot pedal 1, ON/OFF foot pedal with plastic casing (ON by keeping the pedal pressed). Only connects to model: **5040-SX / 5040-SX2**

Art. 290010

Foot pedal 7
Same as foot pedal 1, ON/OFF foot pedal with plastic casing. Only connects to model: **4040**

Art. 290005

Foot pedal 3
Digital foot pedal with rocker function in robust, unbreakable plastic design, black, 197 x 120 x 47 mm, with rubber feet, cable length: 165 cm. Only connects to models: **4030, 4030 S, 5030, 5040, 5040-SX** and **5040-SX2**

Art. 290008

Foot pedal 4
Like foot pedal 3 with rocker function and plastic housing. Only connects to model: **4030-SX, 4030 SX2**

Accessories, replacement filter bags

Replacement micro-fleece filter bags for model 4030 S, 4030-SX, 4030 SX2, 4040 (subject to technical and colour changes)
in a pack of 10 **Art. 195870**
in a pack of 20 **Art. 195871**
in a pack of 50 **Art. 195872**

Accessories for the UVL-36, UVL-36 S

Art. 330904 for UVL-36 S

UV lamp tube 1 pc
(9 W, λ = 365 nm, base = G23)

Art. 330905 for UVL-36 S

UV lamp tube SET OF 4 / 4 pcs.
(9 W, λ = 365 nm, base = G23)

Accessories for nailfan

Art. 230800

Filter unit – washable

Art. 230803

Arm support

Art. 230804

Mesh grid

Accessories for mini nailfan

Art. 295001

Filter (see Nailfan Mini instructions for cleaning)

Accessories for Nailfan Mini-X

Art. 295100

NEW

Filter (see Nailfan Mini-X instructions for cleaning)
2 pcs./pack

Accessories for AC-4000

Art. 334001

Filter set

Art. 334002

UVC lamp

Accessories for TENS/EMS



Art. 359174
Neck electrode
self-adhesive, 1 pc./pack; approx. 13 x 11 cm



Art. 359173
Butterfly electrode
self-adhesive, 2 pcs./pack; 15 x 9.5 cm



Art. 359170
Tissue electrode
self-adhesive, 4 pcs./pack; 40 x 40 mm



Art. 359171
Tissue electrode
self-adhesive, 4 pcs./pack; 45 x 45 mm



Art. 359172
Tissue electrode
self-adhesive, 4 pcs./pack; 45 x 35 mm



Art. 359178
Tissue electrode
self-adhesive, 4 pcs./pack; 45 x 80 mm



Art. 359177
Tissue electrode
self-adhesive, 4 pcs./pack; Ø 32 mm



Art. 359070
Carbon permanent electrode
4 pcs./pack + 4 pcs. self-adhesive gel pads;
40 x 40 mm (reorder of individual gel pads, art. no. 359280)



Art. 359280
Self-adhesive gel pads
for replacing used gel pads for the carbon permanent electrodes (art. no. 359070)
40 x 40 mm, 10 pcs./pack, self-adhesive

Accessories for TENS/EMS



Art. 359910
TENS contact gel
for refreshing the contact surfaces of the electrodes – 50 g tube



Art. 359013
Contact cable
Contact cable for **EMT-4**, in pairs



Art. 359014
Contact cable
Contact cable for **EMT-6**, **T-6** and **IT-6**, in pairs



Art. 359011
Contact cable
Contact cable for **Promed digi 3000**, **Promed tens 1000 s**, in pairs



Art. 359130
Vaginal probe
Special electrode for vaginal use for the treatment of incontinence and for the training of the pelvic floor muscles (women).



Art. 359131
Anal probe
Special electrode for anal use for the treatment of incontinence and for the training of the pelvic floor muscles.



Art. 359990
Battery pack
Battery pack for **Promed tens 1000 s**, **Promed digi tens 3000**, 8.4 V



Art. 359991
9-volt battery pack
9-volt battery pack for **Promed tens 1000 s**, **Promed digi tens 3000**



Art. 359995
Battery, pack of 4
Pack of 4 batteries for Promed **EMT-4**



Art. 359994
Pack of 4 batteries
Battery, pack of 4, for **EMT-4**, **EMT-6**, **T-6**

Therapy accessories


Art. 270960

Swiss Depil – replacement discs
Replacement depilation disc set (6 pcs.)


Art. 402091

Cuff for PBM-3.5, universal (22 - 42 cm)
Arm cuff, universal (22 - 42 cm), grey, for blood pressure monitor


Art. 402093

Power adapter
6-volt mains adapter for blood pressure monitor **PBM-3.5**, **PBW-5.2** and **BDS-700**, **PBM-4.2**

Inhaler accessories


Art. 422001

Protective cap for INH-2.1 and INH-2.1R


Art. 422007

Medicine holder INH-2.1
Spare medicines holder for reordering


Art. 422003

Mouthpiece for INH-2.1 and INH-2.1R
100% latex-free


Art. 422004

Size S mask for INH-2.1 and INH-2.1R
suitable for children, with connecting piece, 100% latex-free


Art. 422005

Size L mask for INH-2.1 and INH-2.1R
suitable for adults, with connecting piece, 100% latex-free


Art. 422006

Adapter for INH-2.1 and INH-2.1R
(image serves only as an example)

Oral hygiene accessories



Art. 470003 Standard nozzle white (white ring)
4043641400873



Art. 470004 Standard nozzle white (grey ring)
4043641400880



Art. 470005 Standard nozzle white (black ring)
4043641400897



Art. 470006 Periodontal nozzle white
4043641400903



Art. 470007 Tongue coating nozzle white
4043641400910



Art. 470008 Orthodontic nozzle white
4043641400927



Art. 470009 Standard nozzle black (black ring)
4043641400934



Art. 470010 Standard nozzle black (grey ring)
4043641400941



Art. 470011 Standard nozzle black (white ring)
4043641400958



Art. 470012 Periodontal nozzle black
4043641400965



Art. 470013 Tongue coating nozzle black
4043641400972



Art. 470014 Orthodontic nozzle black
4043641400989

Accessories for blood glucose monitors


Art. M79039

Control solution for MediTouch 2/MediTouch 2 connect
EAN 40 15588 79039 3


Art. M79038

Test strips for MediTouch 2/MediTouch 2 connect (2 x 25 pcs.)
EAN 40 15588 79038 6


Art. M79002

Lancing device
EAN 40 15588 79002 7


Art. M79028

Lancets for MediTouch 2/MediTouch 2 connect
EAN 40 15588 79028 7

Art. M79049

USB cable for MediTouch 2
EAN 40 15588 79049 2

A

Aroma diffusers:

AL-150 WS aroma diffuser	p. 52
AL-240 FL aroma diffuser	p. 51
AL-400 CL aroma diffuser	p. 51
AL-1300 WS aroma diffuser	p. 50
AL-305 aroma diffuser	p. 52

Aromatic essences/aromatic oils:

Apple-cinnamon (10 ml)	p. 55
Be Happy (10 ml)	p. 54
Eucalyptus (10 ml)	p. 54
Grapefruit (10 ml)	p. 55
Lavender (10 ml)	p. 54
Lemongrass (10 ml)	p. 54
Lime (10 ml)	p. 55
Orange (10 ml)	p. 54
Peach (10 ml)	p. 55
Relax (10 ml)	p. 54
Rose (10 ml)	p. 55
Vanilla (10 ml)	p. 55
Wild cherry (10 ml)	p. 55
Cinnamon (10 ml)	p. 54
Cinnamon orange (10 ml)	p. 55
Lemon (10 ml)	p. 55

B

Blood pressure monitors:

PBM-3.5 upper arm blood pressure monitor	p. 31
HGP-30 wrist blood pressure monitor	p. 31

Blood glucose monitor:

Meditouch 2	p. 37
Meditouch 2 connect DUAL	p. 37

D

Disinfection:

UC-50 ultrasonic cleaner	p. 20
UC-60 ultrasonic cleaner	p. 20

E

Electrical files**Professional and personal care:**

Emotion	p. 18
Feeling	p. 19
Pedimed	p. 16
Pedisenso Duo	p. 17
Sensitive	p. 18
the file - 620	p. 15
the file - 620 deluxe	p. 16
the file - 625	p. 15
the file - 725	p. 12
the file - 825	p. 12
the file - 1030	p. 14
the file - 2020	p. 14
the file - 3030-X	p. 13
the file - 4030-SX2	p. 10
the file - 4040	p. 8
the file - 5040-SX2	p. 11
Ultra pro S	p. 17

I

Inhalation:

INH-2.1 R ultrasonic inhaler	p. 33
------------------------------	-------

Intimate care:

"Orchid" douche	p. 42
"Rose" douche	p. 42

L

Lamps/light:

LTL-749 table lamp	p. 25
LTM-30 table magnifier lamp	p. 25

Air purifiers

AC-4000	p. 27
---------	-------

M

Massage appliances:

ITM Pro hand massage appliance	p. 40
Oral irrigator MD-310	p. 29

P

Pulse oximeters:

PM-200 Pro fingertip pulse oximeter	p. 35
PM-300 OLED pulse oximeter	p. 35
PM-100 connect	p. 36
PM-180	p. 36

S

Grinding bit sets:

Comfort set	p. 76
Callus set	p. 76
Basic modelling set	p. 75
Blue mix modelling set	p. 74
Hard alloy bit set pro modelling set	p. 75
Master modelling set	p. 75
Ceramic modelling set	p. 74
Rainbow modelling set	p. 74
Natural set	p. 76
Personal care set	p. 75

S

Grinding bits:

Blue mix bits	p. 78
Diamond bits	p. 82
Diamond bits with cooling channels	p. 83
Hard alloy bit, blue series	p. 80
Hard alloy bit, yellow series	p. 80
Hard alloy bit, red series	p. 80
Hard alloy bit, specials	p. 81
Hard alloy bit, purple series	p. 79
Ceramic bits	p. 77
Mandrels	p. 85
Polishing bits	p. 86
Rainbow bits	p. 78
Sandpaper bits	p. 85
Sapphire bits	p. 83
Steel milling tool – stainless	p. 81
Stone bits	p. 84
Titanium bits	p. 79

Footwear (bandage shoes):

Flex 2	p. 63
GentleWalk-Hi	p. 67
GentleWalk-Lo	p. 67
Grainau	p. 64
München 2 LXL	p. 68
München 3	p. 68
Pantomobil D	p. 70
Sanicabrio	p. 66
Sanicabrio DXL	p. 65
Sanicabrio LXL	p. 65
Sanisoft	p. 67
Sanisoft D	p. 66
TheraLight 1	p. 70
Theramed D1	p. 69
Theramed D3	p. 69
Wallgau	p. 64

Soles:

Replacement insoles	p. 71
Huber heel spur	p. 73
Huber joint relief	p. 73
Huber hallux-valgus	p. 72

Foot spa:

FB-100 foot spa	p. 41
-----------------	-------

Dust suction:

Nailfan	p. 21
Nailfan mini	p. 21
Nailfan mini- X	p. 22

T

TENS/EMS appliance:

EMT-4	p. 39
EMT-6 combination appliance	p. 38
IT-6	p. 38
Tens 1000 s	p. 39

Thermometers:

IRT-80 infra-red thermometer	p. 34
PFT-3.7 digital thermometer	p. 34

U

UV lamps:

UVL-36 S	p. 24
UVL-54 "All-In"	p. 24
UVL-54-X	p. 23

W

Heating products:

HGP-1.7 heated belt	p. 47
HKP-1.2H cosy heating pad	p. 48
HKP-1.6 XL cosy heating pad	p. 48
KFS-100Foor spa	p. 49
HKP-1.8 XXL heating pad	p. 45
KHP-2.3 cosy heating blanket	p. 44
KHP-2.3G cosy heating blanket	p. 44
KHP-2.3S knitted heating blanket	p. 45
NRP-2.0 shoulder and neck heating pad	p. 46
NRP-2.4 back and neck heating pad	p. 47
NRP-5.4 back and neck heating pad	p. 46
UBS-2.4 underbed	p. 49

Z

Accessories:

Batteries / rechargeable batteries	p. 93
Electrodes	p. 92
Replacement filter bag	p. 90
Foot pedals	p. 90
Handpieces	p. 89
Inhaler	p. 94
Contact gel	p. 93
Contact cable	p. 93
Air purifier accessories	p. 91
Cuffs	p. 94
Accessories for nailfan accessories	p. 91
Accessories for nailfan mini	p. 91
Accessories for nailfan Mini-X	p. 91
Oral irrigator accessories	p. 95
Professional accessories	p. 88
Accessories for the UVL-36, UVL-36 S	p. 91
Vaginal/anal probe	p. 93

 Medisana accessories:

p. 95

We will be happy to send you suitable advertising material:

We continuously update our range of advertising material. Ask our staff team.

Flyer:
(DIN long – different languages)

Electrotherapy flyer (DE)
Art. 993516

Manicure/pedicure flyer (DE)
Art. 991218

Medisana flyer (DE)
Art. M99983

Shoes flyer (DE)
Art. 995114

“Pedibelle” shoes flyer (DE, EN)
DE: Art. 995401; EN: Art. 995402

Heating pads flyer (DE, EN, FR)
DE: Art. 99977; EN: Art. 994501; FR: Art. 994502

Personal care flyer (DE)
Art. 991216

Professional files flyer (DE, EN)
DE: Art. 992020; EN: Art. 992021

Accessories flyer (DE, EN)
DE: Art. 991901; EN: Art. 991902

Blood pressure monitors flyer (DE)
DE: Art. 994101

(Carton) Counter display for DIN long flyers Art. 999981



Your contact persons

Head office / order processing:
Tel.: +49 (0)8821 / 9621-0
info@promed.de

**Management,
marketing and sales:**
Andreas Schuster
Tel.: +49 (0)8821 / 9621-0

Export and mail order retail:
Marsha Lo
Tel.: +49 (0)8821 / 9621-12

Service/returns:
Tel.: +49 (0)8821 / 909 65 50

Field service

Bernd Hegeler
DE region postcode: 2, 3, 49
Tel.: 0151 / 174 815 13

Main office
DE region postcode: 4 (without 49); 5 without 55
Tel.: +49 (0)8821 / 9621-0

Gabriele Zimmer
DE region postcode: 55, 6, 7, 88, 89
Tel.: 0171 / 932 18 85

Ms Gedak
DE region postcode: 80-87; 90-97
Tel.: +49 (0)8821 / 9621-0

Business hours:

Monday - Thursday:
8am - 12pm and 1pm - 5pm
Friday: 8am - 12pm

Returns address

(packages/returns):
Promed GmbH, Bischofstr. 1, D-82490 Farchant

Legal notice

Amendments / GTC

We reserve the right to make technical changes, further developments, printing errors, typesetting errors and mistakes without notice. Our current terms and conditions are available on our website at all times.

Copyright

The information published in the catalogue (text and image material, graphics) and its structure are subject to the copyright of Promed GmbH, and in some cases also to the copyright of third parties. Duplication, dissemination, reproduction and transmission, storage or other use are expressly prohibited without our written consent.

© Copyright 2026, Promed GmbH, Germany



Promed GmbH
Lindenweg 11
D-82490 Farchant
Germany

Tel. +49 (0)8821 / 9621-0
Fax +49 (0)8821 / 9621-21

Email: info@promed.de
Internet: www.promed.de

Managing directors: Andreas Schuster, Thomas Teckentrup
HRB: 135838
Munich Local Court
VAT ID: DE 1283 / 68339

Quick guide (EN):
(DIN A4, in cellophane)

**for TENS/EMS
EMT-4 device**
Art. 983529

**for TENS/EMS
EMT-6 device**
Art. 983530

**Floor display
for Promed products**
Art. 999901

Promed fabric bag
Art. 999986

Ballpoint pens
Art. 940655

A6 writing pad
Art. 999980



www.promed.de



Certified according to DIN EN ISO 13485

Valid from 1 January 2026