

HEALTH & BEAUTY

COMPLETE CATALOGUE EN 1/2026

Special shoes, soles and accessories

WE MAKE YOU MOBILE!

HOLISTIC QUALITY SHOE CONCEPT

Special shoes for short and long-term care are an essential part of the Promed product range.

The continuous further development in collaboration with orthopaedists, medical specialists, rehabilitation centres, geriatric nurses and diabetes specialists has resulted in one of the most holistic collections for the areas of therapy, rehabilitation and post-operative care.

Medical studies document the high standard of the individual models. The best functionality, high-quality materials and first-class workmanship on the one hand, and innovative ideas and pioneering technology on the other, characterise this holistic, high-quality shoe concept.



Softstep sole concept bandage shoes

In most cases, the outer sole of a bandage shoe doesn't look very spectacular. In most cases, this is not the intention anyway, as the focus is on supporting the healing process. This is precisely why, working together with specialists, Promed has come up with a range of intelligent solutions that can provide both relief for patients and support the healing processes.



Front area without tread

This makes everyday walking easier for people suffering from neuropathy and limping.



Metashank inlay

The special Metashank inlay in the forefoot area relieves the pressure on the metatarsal bones and makes every step easier.




Cylindrical cavities

distributed over the entire sole ensure the best cushioning effect. This provides excellent pressure relief not only in the heel, but also in the mid-foot and forefoot areas. It also helps protect the joints and the spine.




Wider sole

Makes everyday walking easier for people suffering from neuropathy and limping.



Diamond pattern tread

together with the "U-I" verge, ensures slip resistance indoors and outdoors, for harmonious heel-to-toe motion.



Bevelled heel area

The heel area, which is bevelled by 20°, ensures a safe tread and ideal heel-to-toe motion.



Ergonomically moulded outer sole (foot shaped)

- Suitable for sanding
- Orthopaedically adjustable
- Made of high-quality PU




THE PROMED SOFTSTEP SOLE CONCEPT COMBINES SAFETY AND COMFORT TO GREAT SUCCESS.

Promed – excellent foot temperature

AIRMESH – THE TEMPERATURE-REGULATING UPPER MATERIAL

[Meshline upper material]



AirMesh with EasyDry lining

The 3-layer shaft material (EasyDry, CosyFit, Mesh) ensures that your feet are comfortable and the temperature is regulated, especially in the summer. The EasyDry fibre, which originates from the medical field, transports

moisture away from the foot and dries quickly. Also particularly suitable for allergy sufferers.

Models: TheraLight 1, TheraLight 2-S



AirMesh + Larina new wool

The 3-layer shaft material (new wool, CosyFit, mesh) ensures warm feet without allowing them to overheat, especially in the winter.

Models: Gentlewalk-Hi, Gentlewalk-Lo




The Promed "AirMesh" is particularly breathable due to its mesh structure and the 3 mm open-pore CosyFit lining foam, which makes it temperature-regulating and ensures a great feeling of well-being for your feet.

It originates from the medical field and is free of toxic components. It is particularly suitable for allergy sufferers and diabetics.



COSYFIT – THE BREATHABLE SHAFT MATERIAL

[Masterline insole]



Models:

- all Masterline, apart from München
- all Classicline, apart from München
- Classic 1 and 2.



This material ensures an exceptionally high level of wear comfort and better pressure relief in sensitive areas.

The multi-layer material absorbs the moisture of the foot well and dries quickly after wearing, thereby ensuring a good foot temperature.



In most cases, corns consist of a painful callus with central, deep-pointing (nail-shaped) spores that are caused by chronic pressure on areas where the skin is near the bone.



Pain in the forefoot can have many causes. The most common cause of such pain is inflammation, often due to incorrect footwear and/or overburdening. It can often occur due to foot deformities that result from accidents.



The various types of inflammation often involve joints, tendons, muscle strands or the periosteum. In the case of a plantar fasciitis, the connective tissue in the sole of the foot is affected.



Hammer toe is a toe deformity in which the toe joints contract when they are bent.



A callus is caused by pressure or friction. This results in a slight inflammation of the skin. To protect the underlying cell layers, the skin is locally stimulated so as to thicken the outer epidermis.



Diabetic foot syndrome, also known colloquially as "diabetic foot", is a syndrome which is associated with diabetes, and which occurs most frequently in patients with type 2 diabetes.



Arthritic foot is the term used to describe an inflammation of the joints in the foot or the toes.



The term splayfoot refers to a spreading of the bones of the metatarsus.



Toe deformities are pathological changes in shape that are either congenital or acquired over the course of life. Increased mechanical pressure can lead to painful calluses and inflammation.



Hallux valgus is the medical term for an X-shaped bend to the big toe in the metatarsal joint. The condition almost always occurs together with splayfoot.



The multi-layer **shaft material** absorbs the moisture of the foot well and dries quickly after wearing, thereby ensuring an **excellent foot temperature**.



Particularly suitable for pressure sensitivity of the feet, especially for **diabetes** care and anyone who wants to support healing processes on their feet.



Particularly suitable for rheumatic symptoms such as rheumatic foot changes, claw and hammer toes, hallux rigidus.



The **Softstep** sole concept, with its cylindrical cavities and neuropathic surface (toe area), combines safety and comfort like never before.



The **AirMesh** shaft material is particularly breathable and temperature-regulating thanks to its mesh structure and open-pore CosyFit foam lining. Also particularly suitable for allergy sufferers and diabetics.



The **"washable insole"** principle allows the sole to be changed at any time and washed at 30 degrees. This prevents unpleasant odours.



Natural new **sheep's wool** ensures comfortably warm feet and a wonderfully soft wearing comfort throughout the shoe.



Upper material made from (perforated) **stretch** for especially good breathability. Adapts to almost any foot shape.



Our shoes are manufactured in the EU under strict quality control



This is a mark of the highest quality. The insoles are the substructure of a shoe. The insole creates a durable, dimensionally stable structure, also without an outer sole, which can therefore be provided with different characteristics.



The hard-wearing material, which is visually similar to suede, ensures a healthy breathable temperature in the shoe. It is abrasion-resistant, tear-resistant and easy care.

General information on the range of Promed shoes:



Machine wash at 30°C



Not suitable for tumble dryer



Do not iron



Do not dry-clean



Do not use bleach

The shoes can be disinfected using any commercially available or common household disinfectant.

SYMBOLS



Ideal for orthopaedic use



Washable



Polyethylene velour



SuperComfort sole wider at the forefoot, protection against tripping, slip resistance



Washable at 30 degrees



Cosyfit lining



Flexible foam insert



Micro velour



Terry lining

Overview of shoes and indications



Designation of shoe	Gentle Walk-Hi	Gentle Walk-Lo	München 2 LXL	München 3	Wallgau	Grainau	Flex 2
Picture of the product							
Art. no.	510341	510331	510230	521241 (b.) 521242 (s.)	520210	520220	520230
HiMi no. (Pair)	31.03.03.4058	31.03.03.4059	31.03.03.4059	31.03.03.4151	31.03.03.4152	31.03.03.4149	31.03.03.4150
Page in the 2026 Promed catalogue	67	67	68	68	64	64	63
Main areas of application							
For general problems with the toes and forefoot				✓	✓	✓	✓✓
Ideal for diabetes care	✓✓	✓✓		✓✓	✓✓	✓✓	✓✓
Ideal for rheumatism sufferers	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓
Ideal for neuropathies	✓	✓	✓	✓	✓	✓	✓
Ideal for wheelchair users	✓✓	✓✓	✓	✓✓	✓✓	✓✓	
For almost any foot shape e.g. hallux, hammer toes, arthritis, swollen feet			✓✓		✓	✓	✓
For inflammations, swellings or bandaged feet			✓		✓	✓	✓✓
For foot deformations			✓✓	✓✓	✓✓	✓✓	✓✓
For use after surgery					✓		
For rheumatic feet	✓✓	✓✓		✓	✓	✓	✓
For feet under heavy strain (e.g. obesity)							
For sore feet e.g. traumatic injuries					✓	✓	✓
For pressure sensitivity, inflammations, swellings and lymphatic feet	✓✓	✓✓			✓		✓
Sores in the forefoot area and on the sole of the foot caused by diabetes							
Sores in the forefoot area and on the sole of the foot caused by other indications such as traumatic injuries							
Sores in the rear-foot area and on the sole of the foot caused by diabetes							
Sores in the rear-foot area and on the sole of the foot caused by other indications such as traumatic injuries							

Bandage shoes ✓ = suitable; ✓✓ = particularly suitable



Theramed D1	Theramed D3	Pantomobil D	TheraLight 1	Sanicabrio DXL	Sanicabrio	Sanicabrio LXL	Sanisoft D	Sanisoft
510160	510180	510140	510311	510130	521120	510150	510170	521110
31.03.03.4042	31.03.03.4041	31.03.03.4043	31.03.03.4059	31.03.03.4045	31.03.03.4020	31.03.03.4046	31.03.03.4047	31.03.03.4022
69	69	70	70	65	66	65	66	67
✓✓	✓✓	✓✓	✓✓	✓✓	✓	✓	✓	✓
					✓✓			✓✓
✓✓	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓
✓	✓	✓	✓	✓	✓✓	✓	✓	✓✓
				✓	✓✓	✓	✓	✓✓
✓	✓✓	✓	✓	✓✓		✓✓		
✓✓	✓✓			✓✓		✓✓	✓✓	
✓✓	✓	✓		✓✓		✓✓	✓✓	✓
✓	✓			✓		✓	✓	
					✓✓			✓✓
✓			✓	✓		✓	✓	
✓	✓	✓		✓✓	✓	✓✓	✓	✓✓
✓✓	✓✓	✓✓		✓	✓	✓✓	✓✓	✓
✓✓	✓✓		✓	✓✓		✓✓	✓✓	
				✓✓	✓✓	✓✓	✓✓	✓✓
	✓✓	✓✓		✓		✓	✓	
		✓✓		✓✓			✓	

SHOE DIMENSIONS PROMED MASTERLINE

Model			SaniCabrio DXL, LXL, SaniSoft D			Theramed D 1			Theramed D 3, Pantomobil D			München 2 LXL		
Design			closed toe / closed heel			open toe / closed heel			open toe / open heel			closed shoe		
English sizes *)	German / French sizes	Foot length in mm	Outsole in mm	Insole in mm	Ball size in mm oversize „J“	Outsole in mm	Insole in mm	Ball size in mm oversize „J“	Outsole in mm	Insole in mm	Ball size in mm oversize „J“	Outsole in mm	Insole in mm	Ball size in mm oversize „J“
3½	36	229	253	238	232	249	239	232	243	239	232	243	234	232
4½	37 ⅓	237	256	242	237	253	243	237	249	244	237	253	238	237
5	38	241	264	249	240	260	250	240	253	248	240	256	242	240
6	39 ⅓	250	273	258	245	264	255	245	260	257	245	264	249	245
6½	40	254	277	261	247	269	259	247	264	262	247	273	258	247
7½	41 ⅓	262	285	268	252	277	268	252	273	271	252	277	261	252
8	42	267	295	276	255	285	275	255	280	278	255	285	268	255
9	43 ⅓	275	298	280	260	295	283	260	285	281	260	295	276	260
9½	44	279	307	290	262	298	288	262	295	289	262	298	280	262
10½	45 ⅓	288	312	294	267	307	296	267	298	295	267	307	290	267
11	46	292	315	300	270	312	301	270	307	303	270	312	294	270
12	47 ⅓	301	324	311	275	318	310	275	312	307	275	315	300	275

SHOE DIMENSIONS PROMED MESHLINE

Model			Gentlewalk Hi SaniCabrio Wallgau, Grainau			Gentlewalk Lo SaniSoft München 3			TheraLight I			Flex 2		
Design			closed toe, closed heel			close toe, close heel			open toe, closed heel			closed toe, open heel		
English sizes *)	German / French sizes	Foot length in mm	Outsole in mm	Insole in mm	Ball width in mm comfort width "G"	Outsole in mm	Insole in mm	Ball width in mm comfort width "G"	Outsole in mm	Insole in mm	Ball width in mm comfort width "G"	Outsole in mm	Insole in mm	Ball width in mm Oversize "J"
3½	36	229	247	238	222	247	238	222	247	235	222	247	238	222
4½	37 ⅓	237	254	242	227	254	242	227	250	240	227	254	242	227
5	38	241	263	249	230	259	246	230	259	247	230	263	249	230
6	39 ⅓	250	268	254	235	273	258	235	263	249	235	268	254	235
6½	40	254	279	261	237	279	261	237	273	258	237	279	261	237
7½	41 ⅓	262	284	268	242	284	268	242	281	270	242	284	268	242
8	42	267	293	276	245	291	275	245	291	278	245	293	276	245
9	43 ⅓	275	299	280	250	299	280	250	291	278	250	299	280	250
9½	44	279	306	290	252	306	290	252	291	278	252	306	290	252
10½	45 ⅓	288	311	294	257	311	294	257	311	300	257	311	294	257
11	46	292	320	303	260	313	300	260	313	302	260	320	303	260
12	47 ⅓	301	321	311	265	321	311	265	321	312	265	321	311	265

*) 1 English stitch = ⅓ inch, i.e. 0.84666667 cm



Art. 520230 (in pairs) | HiMi no.: 31.03.03.4145

Flex 2

Unique bandage shoe shaped as slip-on shoe



- Comfortable slip-on shoe for indoor and outdoor use
- Variable heel strap for 3 possibilities to wear
- Improvable fit due to variable heel and instep strap
- Functioning of the heel strap: pivotable or removable with two press studs, heel width adjustable with Velcro fastener
- Broadened forefoot for comfortable fit and secure footing

- Instep strap with Velcro fastener for perfect and infinite adjustment of the instep height
- Ergonomically shaped, orthopaedically adjustable outer sole with neuropathic zone in the forefoot
- Washable at 30°
- Sizes: 36 - 47, comfortable width J
- Colour: black

Art. 520210 (in pairs) | HiMi no.: 31.03.03.4152

Wallgau Comfortable low bandage shoe Therapeutic comfort shoe

High-quality, 2-layer insole for soft step comfort

- Toe box made of comfortable stretch material
- Particularly suitable for pressure-sensitive feet
- Visually attractive upper material made of elastane (stretch) in the forefoot for particularly good breathability; leather look at rear
- Ergonomically shaped, orthopaedically adjustable outer sole with neuropathy zone in forefoot
- Stretch material in the forefoot area allows for a soft and flexible adjustment to almost any foot shape (hallux, hammer toes, arthritis, swollen feet, etc.)
- Ideal for rheumatic feet and for diabetes care
- Washable at 30 °C
- Sizes: 36 - 47, comfortable width G, colour: black



Art. 520220 (in pairs) | HiMi no.: 31.03.03.4149

Grainau Therapeutic comfort shoe

High-quality, 2-layer insole for soft step comfort

- Toe box made of comfortable stretch material
- Particularly suitable for pressure-sensitive feet
- Visually attractive upper material made of elastane (stretch) in the forefoot for particularly good breathability; leather look at rear
- Ergonomically shaped, orthopaedically adjustable outer sole with neuropathy zone in forefoot
- Stretch material in the forefoot area allows for a soft and flexible adjustment to almost any foot shape (hallux, hammer toes, arthritis, swollen feet, etc.)
- Ideal for rheumatic feet and for diabetes care
- Washable at 30 °C
- Sizes: 36 - 47, comfortable width G, colour: black




Art. 510130 (in pairs) | HiMi no.: 31.03.03.4045

Sanicabrio DXL Over-ankle, extra wide special shoe



- With wide-opening "Cabrio" fastening system at the forefoot and with Velcro fastener at the heel for an easy time putting on the shoe
- Soft and flexible adjustment
- Ideal for swollen feet, hammer toes, for use after operations and with plaster casts and bandaged feet
- Sizes: 36 - 47, extended width J
- Colour: navy blue



Art. 510150 (in pairs) | HiMi no.: 31.03.03.4046

Sanicabrio LXL Low extra-wide special shoe



- With wide-opening "Cabrio" fastening system at the forefoot.
- Easy to put on thanks to Velcro at the heel
- Low shaft at the ankle
- Ideal for solving foot problems such as bandaged feet, heavily stressed, pressure-sensitive feet, hallux and hammer toes
- Sizes: 36 - 47, extended width J
- Colour: navy blue



Wide-opening "Cabrio" fastening system at the forefoot and with Velcro fastener at the heel for an easier time putting on the shoe



Wide-opening "Cabrio" fastening system at the forefoot and with Velcro fastener at the heel for an easy time putting on the shoe



Art. 521120 (in pairs) | HiMi no.: 31.03.03.4020

Sanicabrio

Therapy and rehabilitation shoe with Velcro at the forefoot and heel



- Individual adjustment options
- Also suitable as a first aid shoe
- Ideal for extreme pressure sensitivity and severe swelling
- Dipped edge in the sole for stable hold of soft foam insoles, for example
- Terry lining
- Sizes: 36 - 47, comfortable width G
- Colour: navy blue



Wide-opening "Cabrio" fastening system at the forefoot and with Velcro fastener at the heel for an easy time putting on the shoe



Terry lining



promed

Art. 510170 (in pairs) | HiMi no.: 31.03.03.4047

Sanisoft D

Low-level special shoe



- With large volume and various individual adjustment options in the forefoot and instep area
- A temperature-regulating shoe for diabetics and other foot problems
- Easy to put on and take off, and optimum width adjustment
- Especially suitable for wheelchair users
- Sizes: 36 - 47, extended width J
- Colour: navy blue



Easy to slip on due to wide-opening Velcro fastener in the instep and heel area



promed

Art. 521110 (in pairs) | HiMi no.: 31.03.03.4022

Sanisoft

Therapy and rehabilitation shoe in loafer design



- 2 Velcro fasteners in the forefoot area
- Low shaft cut-out under the ankle
- Dipped edge in the sole for stable hold of soft foam insoles, for example
- Terry lining
- Sizes: 36 - 47, comfortable width G
- Colour: navy blue



Velcro fastener



Terry lining

Lightweight and cushioning sole, slip-resistant and slightly wider in the forefoot

Art. 510341 (in pairs) | HiMi no.: 31.03.03.4.058 (Hi)

Art. 510331 (in pairs) | HiMi no.: 31.03.03.4.059 (Lo)

GentleWalk-Hi / GentleWalk-Lo

Fur boots with Airmesh shaft material and new wool lining



- Breathable "AirMesh" upper material with Larina wool lining for a great foot temperature
- 2 wide-opening Velcro fasteners in the forefoot for easy slip-on and a firm hold
- Pressure-free, ideal fit
- Shock-absorbing "SoftStep" outer sole, suitable for asphalt
- Suitable for diabetics and allergy sufferers



(Hi)



(Lo)

- Machine washable at 30 °C
- Sizes: 36 - 47, comfortable width G
- Colour: black

Art. 510230 (in pairs) | HiMi no.: 31.03.03.4054

München 2 LXL

Stretch micro velour loafer for therapy and rehab

promed





The micro-velour shaft offers a soft, individual hold

Inlaid removable insole

Best finish

Increased forefoot volume and 4-fold stretch material in the forefoot area enable a soft and flexible fit for almost any foot shape

- Increased forefoot volume
- 4-fold stretch material at the forefoot allows for soft and flexible adjustment to almost any foot shape (hallux, hammer toes, arthritis, swollen feet, etc.)
- The stretch material follows your every move
- No painful areas of pressure or friction
- Foot temperature friendly model
- Sizes: 36 - 47, extended width J
- Colour: black



Brown – art. 521241 (in pairs) | HiMi no.: 31.03.03.4151

Black – art. 521242 (in pairs) | HiMi no.: 31.03.03.4151

München 3

Low comfort bandage shoe in a stretch-micro velour combination

promed



Centre and heel part made of micro velour for perfect fit

Instep Velcro fastener for extra convenience when putting the shoe on or taking it off

Slightly damping, anti-slip sole

- Instep Velcro fastener for extra convenience when putting the shoe on or taking it off
- Visually attractive upper material made of perforated stretch Dacron for particularly good breathability
- Ergonomically shaped, orthopaedically adjustable outer sole with neuropathic zone in the forefoot
- Particularly suitable for pressure-sensitive feet
- Stretch material in the forefoot area allows for a soft and flexible adjustment to almost any foot shape (hallux, hammer toes, arthritis, swollen feet, etc.)

- Ideal for rheumatic feet and for diabetes care
- Washable at 30 °C
- Sizes: 36 - 47, comfortable width G
- Colour: brown, black

Art. 510160 (in pairs) | HiMi no.: 31.03.03.4042

Theramed D1

Multifunctional 2-Velcro sandal with closed heel

promed





Breathable and easy-care upper material



Can be completely opened with two Velcro fasteners

- With two Velcro fasteners on the forefoot and instep
- Open toe area, instep strap and closed heel area
- Especially recommended for problems in the toe and forefoot area, such as inflammation, swelling and hallux
- Comfortable in the warmer season due to breathable and easy-care upper material
- Sizes: 36 - 47, oversize J
- Colour: navy blue



Art. 510180 (in pairs) | HiMi no.: 31.03.03.4041

Theramed D3

Multifunctional 2-Velcro sandal with open heel

promed



Three Velcro fasteners for individual adjustment in the forefoot, instep and heel areas



- With three Velcro fasteners for individual adjustment in the forefoot, instep and heel areas. Open heel window
- Recommended for pressure sensitivity in the heel area as well as for inflammation and swelling in the heel and toe area
- Temperature-friendly model
- Comfortable in the warmer season due to breathable upper material
- Sizes: 36 - 47, oversize J
- Colour: navy blue

Art. 510140 (in pairs) | HiMi no.: 31.03.03.4043

Pantomobil D

Lightweight, multifunctional, wide-opening special shoe in slipper design



With wide Velcro fastener in the forefoot for a firm hold and individual fit



Breathable cosy fit upper material (three-layer)

- Recommended for pressure sensitivity, swelling or inflammation in the forefoot and/or heel areas
- Foot temperature-friendly model, particularly pleasant in the warmer season due to breathable cosy-fit upper material (three layers)
- Heel strap included
- Sizes: 36 - 47, oversize J
- Colour: navy blue



promed

Art. 510311 (in pairs)

HiMi no.: 31.03.03.4047

TheraLight 1

Rehabilitation and comfort shoe with open toe area and closed heel cap



Putting the shoe on is extremely easy due to Velcro fastener system

Breathable AirMesh upper material with EasyDry lining



- Breathable "AirMesh" upper material with EasyDry lining for an optimal foot temperature
- wide-opening Velcro fastener to make the shoe easy to put on
- Pressure-free, ideal fit
- Shock-absorbing "SoftStep" outer sole, suitable for asphalt
- Suitable for diabetics and allergy sufferers
- Machine washable at 30°
- Sizes: 36 - 47, comfortable width G
- Colour: black

Art. see below (in pairs)

Replacement soft foam insole (double sole)



Replacement insole (double insole) for the Wallgau and Grainau models
 Sizes: 36 - 47, width J
 Art. 550236

- Soft foam insole 1: Cotton EVA (3.6 mm), black
- Soft foam insole 2: Polyester EVA (4 mm), black

Art. 550235 (in pairs)

Replacement soft foam insole



- High-quality insole with soft step comfort
- Terry-cloth cover
- 4 mm EVA sole
- Suitable for every model in the Promed range
- Sizes: 36 - 47, comfortable width G
- Colour: white

Art. see below (in pairs)

Replacement soft foam insole



Sizes: 36 - 47, width J (in pairs)
 Colour: Black
 Art. 550230 (model: Sanicabrio/DXL/LXL/DS, Sanisoft/D)
 Art. 550231 (model: Theramed D1, TheraLight 1)
 Art. 550232 (model: Theramed D3, Pantomobil D, TheraLight 2-S)

- High-quality insole with soft step comfort
- Terry-cloth cover
- 4 mm EVA sole

Art. 550240 (in pairs)

Replacement fur insole



suitable for hand washing

- Easy care insoles for sensitive feet
- High-quality EVA insole with 100% sheep's wool covering, 4 mm thick with soft tread comfort
- Suitable for the models: Gentlewalk-Hi, Gentlewalk-Lo
- Sizes: 36 - 47, comfortable width G
- Colour: white



HUBER by
Fußtechnologie

promed

THE RIGHT FOUNDATIONS MAKE AL THE DIFFERENCE

Art. 552055

Hallux Valgus

To prevent any further bending of the big toe



HUBER
Fußtechnologie




Excellent for the following symptoms:

- Prevention of hallux valgus
- Sharp pain during the heel-to-toe movement phase of the foot
- Dull pain without stress
- Local pressure sensitivity

Particularly suitable for:

- Sufferers of splayfoot and/or a big toe deformity
- Sufferers of hallux valgus



Features:

- Moisture-absorbing, quick-drying, skin-friendly
- Toxin-free soft microfibre material with pleasant wearing comfort

Sizes (1 pair): 35 - 49

Art. 552050

Heel spur

Heel with U-shaped padded structure




Excellent for the following symptoms:

- Sharp pain when putting weight on the foot
- Dull pain without stress
- Local pressure sensitivity

Particularly suitable for:

- Heel spur, stress on the heel spur
- Persons in occupations that require lots of walking or standing on your feet
- Allergy sufferers and diabetics



Features:

- Moisture-absorbing, quick-drying, skin-friendly
- Toxin-free soft microfibre material with pleasant wearing comfort

Sizes (1 pair): **XS**/35-37, **S**/38-40, **M**/41-43, **L**/44-46, **XL**/47-49



HUBER
Fußtechnologie

Art. 552040 Joint relief sole 1 (up to 80 kg)

Art. 552045 Joint relief sole 2 (over 80 kg)

Joint relief (</> 80 kg)

For targeted, weight-specific pressure relief




Excellent for the following symptoms:


- To relieve stress on the spine and joints
- Complaints regarding the metatarsals (metatarsalgia)
- Burning sole/heel pain, joint problems (ankles, knees, hips)

Particularly suitable for:

- Persons in occupations that require lots of walking or standing on your feet
- Painful, burning, overloaded, tired, pressure-sensitive feet



Helps prevent heavy, tired legs



Helps prevent sweaty feet



Helps prevent burning feet

Features:

- Moisture-absorbing, quick-drying, skin-friendly
- Toxin-free soft microfibre material offering pleasant wearing comfort

Sizes (1 pair): **XS**/35-37, **S**/38-40, **M**/41-43, **L**/44-46, **XL**/47-49



HUBER
Fußtechnologie



www.promed.de



Certified according to DIN EN ISO 13485

Valid from 1 January 2026